

# Product Datasheet

## BAP1 Antibody (BAP1/2665)

### NBP3-07395-100ug

Unit Size: 100 ug

Store at 4C.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-07395](http://www.novusbio.com/NBP3-07395)

Updated 7/16/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-07395](http://www.novusbio.com/reviews/destination/NBP3-07395)



**NBP3-07395-100ug**

BAP1 Antibody (BAP1/2665)

**Product Information**

|                                |                          |
|--------------------------------|--------------------------|
| <b>Unit Size</b>               | 100 ug                   |
| <b>Concentration</b>           | 0.2 mg/ml                |
| <b>Storage</b>                 | Store at 4C.             |
| <b>Clonality</b>               | Monoclonal               |
| <b>Clone</b>                   | BAP1/2665                |
| <b>Preservative</b>            | 0.05% Sodium Azide       |
| <b>Isotype</b>                 | IgG Kappa                |
| <b>Purity</b>                  | Protein A or G purified  |
| <b>Buffer</b>                  | 10 mM PBS with 0.05% BSA |
| <b>Target Molecular Weight</b> | 91 kDa                   |

**Product Description**


|                    |                                                                                                                                                                                                                                                                                              |
|--------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | 200ug/ml of antibody purified from Bioreactor Concentrate by Protein A or G. Prepared in 10 mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0 mg/ml. (NBP3-08506)<br><br>Antibody with azide - store at 2 to 8C. Antibody without azide - store at -20 to -80C. |
| <b>Host</b>        | Mouse                                                                                                                                                                                                                                                                                        |
| <b>Gene ID</b>     | 8314                                                                                                                                                                                                                                                                                         |
| <b>Gene Symbol</b> | BAP1                                                                                                                                                                                                                                                                                         |
| <b>Species</b>     | Human                                                                                                                                                                                                                                                                                        |
| <b>Immunogen</b>   | Recombinant human BAP1 protein fragment (around aa 191-326) (exact sequence is proprietary) (Uniprot: Q92560)                                                                                                                                                                                |

**Product Application Details**

|                              |                                                                                                                |
|------------------------------|----------------------------------------------------------------------------------------------------------------|
| <b>Applications</b>          | ELISA, Protein Array                                                                                           |
| <b>Recommended Dilutions</b> | ELISA, Protein Array                                                                                           |
| <b>Application Notes</b>     | ELISA: For coating, order Ab without BSA.<br>Optimal dilution for a specific application should be determined. |

**Images**

Protein Array: BAP1 Antibody (BAP1/2665) [NBP3-07395] - Analysis of Protein Array containing more than 19,000 full-length human proteins using BAP1 Mouse Monoclonal Antibody (BAP1/2665).





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-07395-100ug**

---

|               |                                                         |
|---------------|---------------------------------------------------------|
| HAF007        | Goat anti-Mouse IgG Secondary Antibody [HRP]            |
| NB720-B       | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP2-57089PEP | BAP1 Recombinant Protein Antigen                        |
| NBL1-07915    | BAP1 Overexpression Lysate                              |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-07395](http://www.novusbio.com/reviews/submit/NBP3-07395)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

