

Product Datasheet

Human PACAP-38 ELISA Kit (Colorimetric) NBP3-08165

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-08165

Updated 10/7/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-08165




NBP3-08165**Human PACAP-38 ELISA Kit (Colorimetric)**

Product Information	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Storage of components varies. See protocol for specific instructions.
Product Description	
Gene ID	116
Gene Symbol	ADCYAP1
Species	Human
Kit Components	Micro ELISA Plate (Dismountable) - 1 EA, Reference Standard - 2 vials, Concentrated Biotinylated Detection Ab (100X) - 120 uL, Concentrated HRP Conjugate (100X) - 120 uL, Reference Standard & Sample Diluent - 20 mL, Biotinylated Detection Ab Diluent - 14
Standard Curve Range	0.31 - 20 ng/mL
Sensitivity	0.19 ng/mL
Inter Assay	%CV < 4.95
Intra Assay	%CV < 6.86
Assay Type	Sandwich-ELISA
Spiking Recovery	85 - 106%
Suitable Sample Type	Serum, Plasma, Cell supernatant
Product Application Details	
Applications	ELISA
Recommended Dilutions	ELISA




Images

ELISA: Human PACAP-38 ELISA Kit (Colorimetric) [NBP3-08165] - Samples were spiked with high concentrations of Human PACAP-38 and diluted with Reference Standard & Sample Diluent to produce samples with values within the range of the assay.



ELISA: Human PACAP-38 ELISA Kit (Colorimetric) [NBP3-08165] - Standard Curve Reference





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-08165

M6000B-1	IL-6 [HRP]
NBP3-28393-100ug	PACAP-38 Antibody (Lilly patent anti-PACAP) - Humanized
NB300-109	Tyrosine Hydroxylase Antibody
MAB1417	Insulin Antibody (182410) [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-08165

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

