

Product Datasheet

p38 gamma/SAPK3 Antibody (R13-L) NBP3-08207

Unit Size: 50 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-08207

Updated 11/12/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-08207




NBP3-08207

p38 gamma/SAPK3 Antibody (R13-L)


Product Information	
Unit Size	50 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	R13-L
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	20 mM Tris-HCl, pH 8.0, 10 mg/ml BSA
Product Description	
Description	This antibody is immunoaffinity purified with immunogenic peptide as a ligand.
Host	Rabbit
Gene ID	6300
Gene Symbol	MAPK12
Species	Human, Mouse, Rat
Immunogen	Peptide derived from the C-terminal sequence of human p38 gamma/SAPK3. Antibody recognizes the epitope located between Arg355 - Pro366.
Notes	This antibody is immunoaffinity purified with immunogenic peptide as a ligand.
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/Immunofluorescence, Immunoprecipitation
Recommended Dilutions	Western Blot 1:1000, ELISA 1:10000-1:20000, Immunocytochemistry/Immunofluorescence, Immunoprecipitation

Images

Western Blot: p38 gamma/SAPK3 Antibody (R13-L) [NBP3-08207] - Western blot analysis of p38 gamma/SAPK3. Lane A - human recombinant 38-alpha(260-360a.a., N-terminal GST tag; H00001432-Q01), lane B - human recombinant p38-beta(with C-terminal DDK tag, NBL1-12870), lane C - human recombinant p38 gamma/SAPK3 (1-368a.a., N-terminal GST tag, H00006300-P01), lane D - human recombinant p38-delta (1-366a.a., N-terminal GST tag, H00005603-P01). 500ng of recombinant protein was loaded in each well.



Immunocytochemistry/Immunofluorescence: p38 gamma/SAPK3 Antibody (R13-L) [NBP3-08207] - Representative picture of p38 gamma/SAPK3 expression in HEK293 cells, visualized with clonal rabbit anti-p38 gamma/SAPK3 monospecific antibody. Primary antibody dilution - 1:300.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-08207

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-72551-100ug	Recombinant Human p38 gamma/SAPK3 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-08207

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

