

# Product Datasheet

## Resistin Antibody (RETN/4326) [DyLight 488] NBP3-08258G

Unit Size: 100 ul

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-08258G](http://www.novusbio.com/NBP3-08258G)

Updated 10/26/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-08258G](http://www.novusbio.com/reviews/destination/NBP3-08258G)



**NBP3-08258G**

Resistin Antibody (RETN/4326) [DyLight 488]

| <b>Product Information</b> |   |
|----------------------------|---|
| <b>Unit Size</b>           | 100 ul  |
| <b>Concentration</b>       | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>             | Store at 4C in the dark.  |
| <b>Clonality</b>           | Monoclonal  |
| <b>Clone</b>               | RETN/4326   |
| <b>Preservative</b>        | 0.05% Sodium Azide  |
| <b>Isotype</b>             | IgG2b Kappa   |
| <b>Conjugate</b>           | DyLight 488   |
| <b>Purity</b>              | Protein A or G purified   |
| <b>Buffer</b>              | 50mM Sodium Borate  |

| <b>Product Description</b> |   |
|----------------------------|---|
| <b>Host</b>                | Mouse   |
| <b>Gene ID</b>             | 56729   |
| <b>Gene Symbol</b>         | RETN  |
| <b>Species</b>             | Human   |
| <b>Immunogen</b>           | Recombinant fragment of human Resistin protein (around aa 19-108) (exact sequence is proprietary) (Uniprot: Q9HD89) |
| <b>Notes</b>               | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                   |

| <b>Product Application Details</b> |  |
|------------------------------------|--|
| <b>Applications</b>                | Western Blot, Immunohistochemistry-Paraffin, Protein Array             |
| <b>Recommended Dilutions</b>       | Western Blot, Immunohistochemistry-Paraffin, Protein Array             |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined. |





### Novus Biologicals USA

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### Bio-Techne Canada

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### Bio-Techne Ltd

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### General Contact Information

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### Products Related to NBP3-08258G

---

|                |  |
|----------------|--|
| NBP1-43317G    | Mouse IgG2b Kappa Light Chain Isotype Control (MG2b) [DyLight 488] |
| NBP2-61343-5ug | Recombinant Rat Resistin Protein                                   |
| 210-TA-005     | TNF-alpha [Unconjugated]   |
| DRSN00         | Resistin [HRP]   |

---

### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-08258G](http://www.novusbio.com/reviews/submit/NBP3-08258G)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

