

Product Datasheet

ZNF408 Antibody (PCRP-ZNF408-1E5) - Azide and BSA Free **NBP3-08324**

Unit Size: 100 ug

Store at -20 to -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-08324

Updated 7/16/2024 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-08324



NBP3-08324

ZNF408 Antibody (PCRP-ZNF408-1E5) - Azide and BSA Free

Product Information	
Unit Size	100 ug
Concentration	1 mg/ml
Storage	Store at -20 to -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	PCRP-ZNF408-1E5
Preservative	No Preservative
Isotype	IgG2b
Purity	Protein A or G purified
Buffer	10 mM PBS
Target Molecular Weight	118 kDa


Product Description	
Description	1.0 mg/ml of antibody purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS WITHOUT BSA & azide. Also available at 200 ug/ml WITH BSA & azide (NBP3-07366). Antibody with azide - store at 2 to 8C. Antibody without azide - store at -20 to -80 C.
Host	Mouse
Gene ID	79797
Gene Symbol	ZNF408
Species	Human
Immunogen	Recombinant full-length human ZNF408 protein (Uniprot: Q9H9D4)

Product Application Details	
Applications	Western Blot, Flow Cytometry, Protein Array
Recommended Dilutions	Western Blot 1-2 ug/ml, Flow Cytometry 1-2 ug/million cells, Protein Array
Application Notes	Optimal dilution for a specific application should be determined.




Images

Flow Cytometry: ZNF408 Antibody (PCRP-ZNF408-1E5) [NBP3-08324] - Flow Cytometric Analysis of PFA-fixed HeLa cells. ZNF408 Mouse Monoclonal Antibody (PCRP-ZNF408-1E5) followed by goat anti- Mouse IgG-CF488 (blue); unstained cells (red).



Protein Array: ZNF408 Antibody (PCRP-ZNF408-1E5) [NBP3-08324] - Analysis of Protein Array containing more than 19,000 full-length human proteins using ZNF408 Mouse Monoclonal Antibody (PCRP-ZNF408-1E5).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-08324

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP2-27231	Mouse IgG2b Isotype Control (MPC-11)
NBP1-82280PEP	ZNF408 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-08324

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

