

# Product Datasheet

## Prohibitin Antibody (PHB/3227) [DyLight 755] NBP3-08336IR

Unit Size: 100 ul

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-08336IR](http://www.novusbio.com/NBP3-08336IR)

Updated 10/26/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-08336IR](http://www.novusbio.com/reviews/destination/NBP3-08336IR)



**NBP3-08336IR**

Prohibitin Antibody (PHB/3227) [DyLight 755]

| Product Information |   |
|---------------------|---|
| Unit Size           | 100 ul  |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.  |
| Clonality           | Monoclonal  |
| Clone               | PHB/3227  |
| Preservative        | 0.05% Sodium Azide  |
| Isotype             | IgG2a Kappa   |
| Conjugate           | DyLight 755   |
| Purity              | Protein A or G purified   |
| Buffer              | 50mM Sodium Borate  |

| Product Description |   |
|---------------------|---|
| Host                | Mouse   |
| Gene ID             | 5245  |
| Gene Symbol         | PHB1  |
| Species             | Human   |
| Marker              | Mitochondrial Marker  |
| Immunogen           | Recombinant human Prohibitin protein fragment (aa167-261) (exact sequence is proprietary) (Uniprot: P35232) |
| Notes               | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                           |

| Product Application Details |  |
|-----------------------------|--|
| Applications                | Western Blot, Immunohistochemistry-Paraffin, Protein Array             |
| Recommended Dilutions       | Western Blot, Immunohistochemistry-Paraffin, Protein Array             |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-08336IR**

---

|              |  |
|--------------|--|
| NBP1-96981IR | Mouse IgG2a Kappa Isotype Control (M2AK) [DyLight 755] |
| NBP2-23364   | Recombinant Human Prohibitin His Protein               |
| 210-TA-005   | TNF-alpha [Unconjugated]                               |
| NBL1-14348   | Prohibitin Overexpression Lysate                       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-08336IR](http://www.novusbio.com/reviews/submit/NBP3-08336IR)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

