

Product Datasheet

GATA-3 Antibody (GATA3/2441) [PE] NBP3-08654PE

Unit Size: 100 ul

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-08654PE

Updated 10/26/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-08654PE



NBP3-08654PE

GATA-3 Antibody (GATA3/2441) [PE]

Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	GATA3/2441
Preservative	0.05% Sodium Azide
Isotype	IgG2b Kappa
Conjugate	PE
Purity	Protein A or G purified
Buffer	PBS

Product Description	
Host	Mouse
Gene ID	2625
Gene Symbol	GATA3
Species	Human
Marker	Breast and Urothelial Marker
Specificity/Sensitivity	GATA-3 is a zinc finger transcription factor and plays an important role in promoting and directing cell proliferation, development, and differentiation in many tissues and cell types. GATA-3 expression is primarily seen in breast and urothelial carcinomas. Therefore, GATA3 antibody can be useful in the identification of unknown primary carcinoma when carcinomas of the breast or bladder are a possibility.
Immunogen	Recombinant fragment of human GATA-3 protein (around aa 357-436) (exact sequence is proprietary) (Uniprot: P23771)

Product Application Details	
Applications	Western Blot, ELISA, Flow Cytometry, Protein Array
Recommended Dilutions	Western Blot, Flow Cytometry, ELISA, Protein Array
Application Notes	Optimal dilution of this antibody should be experimentally determined.

Images

GATA-3 Antibody (GATA3/2441) [PE] [NBP3-08654PE] - Vial of PE conjugated antibody. PE has two excitation maxima, 498 nm excited by the Blue laser (488 nm) and 565 nm excited by the Yellow-Green laser (561 nm). Both result in emission at 578 nm.



PE

LASER (nm)	FILTER
Blue (488)	582/15
Y-G (561)	

EXCITATION MAX (nm)	EMISSION MAX (nm)
498	578
565	



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-08654PE

NBP1-43317PE	Mouse IgG2b Kappa Light Chain Isotype Control (MG2b) [PE]
H00002625-P02-2ug	Recombinant Human GATA-3 GST (N-Term) Protein
210-TA-005	TNF-alpha [Unconjugated]
H00002625-Q02-10ug	Recombinant Human GATA-3 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-08654PE

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

