

# Product Datasheet

## Spastin Antibody (Sp 3G11-1) [Allophycocyanin] NBP3-08681APC

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-08681APC](http://www.novusbio.com/NBP3-08681APC)

Updated 11/11/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-08681APC](http://www.novusbio.com/reviews/destination/NBP3-08681APC)



**NBP3-08681APC**

Spastin Antibody (Sp 3G11-1) [Allophycocyanin]

| Product Information         |                                                                                                                  |
|-----------------------------|------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                           |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                      |
| Storage                     | Store at 4C in the dark.                                                                                         |
| Clonality                   | Monoclonal                                                                                                       |
| Clone                       | Sp 3G11-1                                                                                                        |
| Preservative                | 0.05% Sodium Azide                                                                                               |
| Isotype                     | IgG2a Kappa                                                                                                      |
| Conjugate                   | Allophycocyanin                                                                                                  |
| Purity                      | Protein A or G purified                                                                                          |
| Buffer                      | PBS                                                                                                              |
| Product Description         |                                                                                                                  |
| Host                        | Mouse                                                                                                            |
| Gene ID                     | 6683                                                                                                             |
| Gene Symbol                 | SPAST                                                                                                            |
| Species                     | Human, Mouse, Rat                                                                                                |
| Immunogen                   | Recombinant full-length human Spastin protein. (Uniprot: Q9UP0)                                                  |
| Product Application Details |                                                                                                                  |
| Applications                | Western Blot, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin |
| Recommended Dilutions       | Western Blot, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                           |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to NBP3-08681APC**

|                    |                                                            |
|--------------------|------------------------------------------------------------|
| NBP1-96981APC      | Mouse IgG2a Kappa Isotype Control (M2AK) [Allophycocyanin] |
| H00006683-Q01-10ug | Recombinant Human Spastin GST (N-Term) Protein             |
| NBP1-90230PEP      | Spastin Recombinant Protein Antigen                        |
| NBP2-13075         | beta Amyloid Antibody (MOAB-2) - BSA Free                  |

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-08681APC](http://www.novusbio.com/reviews/submit/NBP3-08681APC)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

