

Product Datasheet

bcl-x Antibody (BCL2L1/2406) - Azide and BSA Free NBP3-08777

Unit Size: 100 ug

Store at -20 to -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-08777

Updated 7/16/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-08777



NBP3-08777


bcl-x Antibody (BCL2L1/2406) - Azide and BSA Free

| Product Information | |
|-----------------------------|--|
| Unit Size | 100 ug |
| Concentration | 1 mg/ml |
| Storage | Store at -20 to -80C. Avoid freeze-thaw cycles. |
| Clonality | Monoclonal |
| Clone | BCL2L1/2406 |
| Preservative | No Preservative |
| Isotype | IgG |
| Purity | Protein A or G purified |
| Buffer | 10 mM PBS |
| Target Molecular Weight | 27 kDa |
| Product Description | |
| Description | <p>1.0 mg/ml of antibody purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS WITHOUT BSA & azide. Also available at 200 ug/ml WITH BSA & azide (NBP3-07455).</p> <p>Antibody with azide - store at 2 to 8C. Antibody without azide - store at -20 to -80 C.</p> |
| Host | Mouse |
| Gene ID | 598 |
| Gene Symbol | BCL2L1 |
| Species | Human |
| Marker | Apoptosis Marker |
| Specificity/Sensitivity | Recognizes a protein of 27kDa, identified as the Bcl-X protein. This monoclonal antibody shows no cross-reaction with Bcl-2 or Bax protein. Bcl-X has two isoforms, Bcl-XL (long), a 241-amino acid protein which suppresses cell death. And Bcl-XS (short), a 178-amino acid protein lacking a 63-amino acid domain which functions as a dominant inhibitor of Bcl-2. |
| Immunogen | Recombinant human bcl-x protein fragment (around aa 1-135) (exact sequence is proprietary) (Uniprot: Q07817) |
| Product Application Details | |
| Applications | ELISA, Immunohistochemistry-Paraffin |
| Recommended Dilutions | ELISA 1-5 ug/ml, Immunohistochemistry-Paraffin |
| Application Notes | ELISA: 1-5ug/ml. Order antibody without BSA. Optimal dilution for a specific application should be determined. |




Images

Immunohistochemistry-Paraffin: bcl-x Antibody (BCL2L1/2406) - Azide and BSA Free [NBP3-08777] - Formalin-fixed, paraffin-embedded Hodgkin's lymphoma stained with bcl-x Mouse Monoclonal Antibody (BCL2L1/2406).



Immunohistochemistry-Paraffin: bcl-x Antibody (BCL2L1/2406) - Azide and BSA Free [NBP3-08777] - Formalin-fixed, paraffin-embedded Hodgkin's lymphoma stained with bcl-x Mouse Monoclonal Antibody (BCL2L1/2406).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-08777

| | |
|----------------|---|
| HAF007 | Goat anti-Mouse IgG Secondary Antibody [HRP] |
| NB720-B | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-97019-5mg | Mouse IgG Isotype Control |
| 210-TA-005 | TNF-alpha [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-08777

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

