

Product Datasheet

DREF Antibody (PCRP-ZBED1-1E1) - Azide and BSA Free NBP3-08861

Unit Size: 100 ug

Store at -20 to -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-08861

Updated 7/16/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-08861



NBP3-08861

DREF Antibody (PCRP-ZBED1-1E1) - Azide and BSA Free

Product Information	
Unit Size	100 ug
Concentration	1 mg/ml
Storage	Store at -20 to -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	PCRP-ZBED1-1E1
Preservative	No Preservative
Isotype	IgG2a
Purity	Protein A or G purified
Buffer	10 mM PBS
Target Molecular Weight	86 kDa


Product Description	
Description	1.0 mg/ml of antibody purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS WITHOUT BSA & azide. Also available at 200 ug/ml WITH BSA & azide (NBP3-07536). Antibody with azide - store at 2 to 8C. Antibody without azide - store at -20 to -80 C.
Host	Mouse
Gene ID	9189
Gene Symbol	ZBED1
Species	Human
Immunogen	Recombinant full-length human DREF protein (Uniprot: O96006)

Product Application Details	
Applications	Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Protein Array
Recommended Dilutions	Western Blot 1-2 ug/ml, Flow Cytometry 1-2 ug/million cells, Immunocytochemistry/ Immunofluorescence 1-2 ug/ml, Protein Array
Application Notes	Optimal dilution for a specific application should be determined.




Images


Immunocytochemistry/Immunofluorescence: DREF Antibody (PCRP-ZBED1-1E1) [NBP3-08861] - Immunofluorescence Analysis of PFA-fixed MCF-7 cells using DREF Mouse Monoclonal Antibody (PCRP-ZBED1-1E1) followed by goat anti- Mouse IgG-CF488 (green). CF640A phalloidin (red).




Flow Cytometry: DREF Antibody (PCRP-ZBED1-1E1) [NBP3-08861] - Flow Cytometric Analysis of PFA-fixed HeLa cells. DREF Mouse Monoclonal Antibody (PCRP-ZBED1-1E1) followed by goat anti- Mouse IgG-CF488 (blue); unstained cells (red).



Immunocytochemistry/Immunofluorescence: DREF Antibody (PCRP-ZBED1-1E1) [NBP3-08861] - Immunofluorescence Analysis of PFA-fixed HeLa cells using DREF Mouse Monoclonal Antibody (PCRP-ZBED1-1E1) followed by goat anti- Mouse IgG-CF488 (green). CF640A phalloidin (red).



Protein Array: DREF Antibody (PCRP-ZBED1-1E1) [NBP3-08861] - Analysis of Protein Array containing more than 19,000 full-length human proteins using DREF Mouse Monoclonal Antibody (PCRP-ZBED1-1E1).



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-08861

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96778	Mouse IgG2a Isotype Control (M2A)
H00009189-P01-10ug	Recombinant Human DREF GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-08861

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

