

# Product Datasheet

## Moesin Antibody (rMSN/492) [PerCP] NBP3-08923PCP

Unit Size: 100 ul

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-08923PCP](http://www.novusbio.com/NBP3-08923PCP)

Updated 10/26/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-08923PCP](http://www.novusbio.com/reviews/destination/NBP3-08923PCP)



**NBP3-08923PCP**

Moesin Antibody (rMSN/492) [PerCP]

| Product Information         |                                                                                                                           |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ul                                                                                                                    |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                               |
| Storage                     | Store at 4C in the dark.                                                                                                  |
| Clonality                   | Monoclonal                                                                                                                |
| Clone                       | rMSN/492                                                                                                                  |
| Preservative                | 0.05% Sodium Azide                                                                                                        |
| Isotype                     | IgG1 Kappa                                                                                                                |
| Conjugate                   | PerCP                                                                                                                     |
| Purity                      | Protein A or G purified                                                                                                   |
| Buffer                      | PBS                                                                                                                       |
| Product Description         |                                                                                                                           |
| Host                        | Mouse                                                                                                                     |
| Gene ID                     | 4478                                                                                                                      |
| Gene Symbol                 | MSN                                                                                                                       |
| Species                     | Human, Rat                                                                                                                |
| Immunogen                   | Recombinant full-length human Moesin protein (Uniprot: P26038)                                                            |
| Product Application Details |                                                                                                                           |
| Applications                | Western Blot, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin |
| Recommended Dilutions       | Western Blot, Flow Cytometry, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                    |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to NBP3-08923PCP**

|                     |                                                       |
|---------------------|-------------------------------------------------------|
| NBP1-43319PCP-0.5ml | Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1) [PerCP] |
| H00004478-P01-10ug  | Recombinant Human Moesin GST (N-Term) Protein         |
| 7268-CT-100         | CTLA-4 [Unconjugated]                                 |
| H00004478-Q01-10ug  | Recombinant Human Moesin GST (N-Term) Protein         |

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-08923PCP](http://www.novusbio.com/reviews/submit/NBP3-08923PCP)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

