

Product Datasheet

TWEAK/TNFSF12 Antibody (CARL-1) - BSA Free NBP3-09144

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-09144

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-09144



NBP3-09144

TWEAK/TNFSF12 Antibody (CARL-1) - BSA Free

Product Information

Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	CARL-1
Preservative	15mM Sodium Azide
Isotype	IgG3 Kappa
Purity	Protein A purified
Buffer	Tris buffered saline (TBS), pH 8.0

Product Description


Description	Novus Biologicals Mouse TWEAK/TNFSF12 Antibody (CARL-1) - BSA Free (NBP3-09144) is a monoclonal antibody validated for use in IHC and Flow. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	8742
Gene Symbol	TNFSF12
Species	Human
Specificity/Sensitivity	The mouse monoclonal antibody CARL-1 recognizes an extracellular epitope of CD255 / TWEAK, a type II transmembrane protein of the TNF superfamily able to induce apoptosis weakly in many cell types.
Immunogen	human TWEAK/TNFSF12-transfected 2PK-3 cells

Product Application Details

Applications	Flow Cytometry, Immunohistochemistry
Recommended Dilutions	Flow Cytometry 1-4 ug/ml, Immunohistochemistry

Images

Flow Cytometry: TWEAK/TNFSF12 Antibody (CARL-1) [NBP3-09144] - Separation of CD255 transfected L5178Y cells (red-filled) from nontransfected L5178Y cells (black-dashed) in flow cytometry analysis (surface staining) stained using anti-human CD255 (CARL-1) purified antibody (concentration in sample 9 ug/ml) GAM APC.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-09144

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96978	Mouse IgG3 Kappa Light Chain Isotype Control (MG3K)
NBP2-46811PEP	TWEAK/TNFSF12 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-09144

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

