

# Product Datasheet

## OR3A3 Antibody - BSA Free NBP3-09739-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-09739](http://www.novusbio.com/NBP3-09739)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-09739](http://www.novusbio.com/reviews/destination/NBP3-09739)



**NBP3-09739-100ul**

OR3A3 Antibody - BSA Free


| Product Information     |  |
|-------------------------|--|
| Unit Size               | 100 ul   |
| Concentration           | 0.5 mg/ml  |
| Storage                 | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| Clonality               | Polyclonal   |
| Preservative            | 0.09% Sodium Azide   |
| Purity                  | Affinity purified  |
| Buffer                  | PBS, 2% Sucrose  |
| Target Molecular Weight | 35 kDa   |

| Product Description |   |
|---------------------|---|
| Description         | Novus Biologicals Rabbit OR3A3 Antibody - BSA Free (NBP3-09739) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host                | Rabbit  |
| Gene ID             | 8392  |
| Gene Symbol         | OR3A3   |
| Species             | Human   |
| Immunogen           | The immunogen is a synthetic peptide directed towards the C-terminal region of Human OR3A3 (NP_036505). Peptide sequence<br>MRLGSVESSDKDKGVGVFMTVINPMLNPLIYSLRNTDVQGALCQLLVGKR        |

| Product Application Details |                        |
|-----------------------------|------------------------|
| Applications                | Western Blot           |
| Recommended Dilutions       | Western Blot 1.0 ug/ml |

**Images**

Western Blot: OR3A3 Antibody [NBP3-09739] - Western blot analysis of OR3A3 in Human 293T Whole Cell. Antibody dilution at 1ug/ml





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-09739-100ul**

---

|                    |   |
|--------------------|---|
| NBP2-33376H        | Blue Marker Antibody (6F4-F6) [HRP]                 |
| HAF008             | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160             | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| H00008392-P01-10ug | Recombinant Human OR3A3 GST (N-Term) Protein        |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-09739](http://www.novusbio.com/reviews/submit/NBP3-09739)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

