

Product Datasheet

EIF3G Antibody - BSA Free

NBP3-10537-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-10537

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-10537



NBP3-10537-100ul

EIF3G Antibody - BSA Free

Product Information	
Unit Size	100 ul
Concentration	0.5 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Purity	Affinity purified
Buffer	PBS, 2% Sucrose
Target Molecular Weight	35 kDa


Product Description	
Description	Novus Biologicals Rabbit EIF3G Antibody - BSA Free (NBP3-10537) is a polyclonal antibody validated for use in IHC, WB and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	8666
Gene Symbol	EIF3G
Species	Human, Mouse
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human EIF3G (NP_003746). Peptide sequence LRDGASRRGESMQPNRRADDNATIRVTNLSEDTRETDLQELFRPFGSISR

Product Application Details	
Applications	Western Blot, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot, Immunohistochemistry, Immunoprecipitation



Images


Western Blot: EIF3G Antibody [NBP3-10537] - Host: Rabbit. Target Name: EIF3G. Sample Type: MCF7. Antibody Dilution: 1.0ug/ml EIF3G is supported by BioGPS gene expression data to be expressed in MCF7



Western Blot: EIF3G Antibody [NBP3-10537] - WB Suggested Anti-EIF3G Antibody Titration: 0.2-1 ug/ml. ELISA Titer: 1:62500. Positive Control: Jurkat cell lysate



Western Blot: EIF3G Antibody [NBP3-10537] - Host: Rabbit. Target Name: EIF3G. Sample Type: 293T. Antibody Dilution: 1.0ug/ml EIF3G is supported by BioGPS gene expression data to be expressed in HEK293T





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-10537-100ul

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
H00008666-Q01-10ug	Recombinant Human EIF3G GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-10537

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

