

Product Datasheet

B7-1/CD80 Antibody (RM80) [mFluor Violet 450 SE] NBP3-11977MFV450

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-11977MFV450

Updated 9/20/2023 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-11977MFV450



NBP3-11977MFV450

B7-1/CD80 Antibody (RM80) [mFluor Violet 450 SE]


Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	RM80
Preservative	0.05% Sodium Azide
Isotype	IgG2a Kappa
Conjugate	mFluor Violet 450 SE
Purity	Protein A purified
Buffer	50mM Sodium Borate

Product Description	
Host	Rat
Gene ID	941
Gene Symbol	CD80
Species	Mouse
Specificity/Sensitivity	This antibody is specific for murine B7-1/CD80, a ~55 kDa member of the Ig superfamily. B7-1/CD80 is expressed by macrophages, dendritic cells, activated B cells, and activated T cells. B7-1/CD80 has high affinity for binding to two T cell surface antigens, CD28 and CD152 (CTLA-4). The interaction of CD28 and CD152 with B7-1/CD80 is crucial in T-B cell communication leading to the activation of T and B cells, respectively.
Immunogen	This antibody was raised by immunising Sprague Dawley rats with a mouse B cell line, BCL1 and fusing immune splenocytes with P3U1 myeloma cells.

Product Application Details	
Applications	Western Blot, Flow Cytometry, Immunohistochemistry, Block/Neutralize, Functional Assay
Recommended Dilutions	Western Blot, Flow Cytometry, Immunohistochemistry, Block/Neutralize, Functional Assay
Application Notes	Optimal dilution of this antibody should be experimentally determined.

Images

B7-1/CD80 Antibody (RM80) [mFluor Violet 450 SE] [NBP3-11977MFV450] - Vial of mFluor Violet 450 conjugated antibody. mFluor Violet 450 is optimally excited at 406 nm by the Violet laser (405 nm) and has an emission maximum of 445 nm.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-11977MFV450

NBP2-26578	Recombinant Mouse B7-1/CD80 Protein
210-TA-005	TNF-alpha [Unconjugated]
740-B1-100	B7-1/CD80 [Unconjugated]
6507-IL-010/CF	IL-4 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-11977MFV450

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications



