

Product Datasheet

4-1BB Ligand/TNFSF9 Antibody (AT113-2) [mFluor Violet 450 SE] NBP3-11987MFV450

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-11987MFV450

Updated 9/20/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-11987MFV450



NBP3-11987MFV450

4-1BB Ligand/TNFSF9 Antibody (AT113-2) [mFluor Violet 450 SE]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	AT113-2
Preservative	0.05% Sodium Azide
Isotype	IgG1 Kappa
Conjugate	mFluor Violet 450 SE
Purity	Protein A purified
Buffer	50mM Sodium Borate

Product Description	
Host	Rat
Gene ID	8744
Gene Symbol	TNFSF9
Species	Mouse
Specificity/Sensitivity	This antibody is specific for 4-1BB ligand (also known as CD137 ligand), which is a member of the TNF superfamily and is a type II membrane protein. 4-1BB ligand is a co-stimulatory molecule that is expressed on antigen presenting cells (DCs, monocytes/macrophages, B cells) and is upregulated upon activation of these cells. Both of the 4-1BB ligand and its receptor, TNFRSF9/4-1BB, are involved in the antigen presentation process and in the generation of cytotoxic T cells. In direct ELISAs, no cross-reactivity with recombinant human RELT is observed.
Immunogen	This antibody was raised by immunising LOU rats with soluble recombinant 4-1BB Ligand/TNFSF9:Fc fusion protein.

Product Application Details	
Applications	Western Blot, Flow Cytometry, Block/Neutralize
Recommended Dilutions	Western Blot, Flow Cytometry, Block/Neutralize
Application Notes	Optimal dilution of this antibody should be experimentally determined.

Images

4-1BB Ligand/TNFSF9 Antibody (AT113-2) [mFluor Violet 450 SE] [NBP3-11987MFV450] - Vial of mFluor Violet 450 conjugated antibody. mFluor Violet 450 is optimally excited at 406 nm by the Violet laser (405 nm) and has an emission maximum of 445 nm.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-11987MFV450

NBP2-26580	Recombinant Mouse 4-1BB Ligand/TNFSF9 Protein
210-TA-005	TNF-alpha [Unconjugated]
1246-4L-025/CF	4-1BB Ligand/TNFSF9 [Unconjugated]
285-IF-100	IFN-gamma [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-11987MFV450

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications



