

Product Datasheet

CD30 Ligand/TNFSF8 Antibody (RM153) [mFluor Violet 450 SE]

NBP3-11991MFV450

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-11991MFV450

Updated 9/20/2023 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-11991MFV450



NBP3-11991MFV450

CD30 Ligand/TNFSF8 Antibody (RM153) [mFluor Violet 450 SE]

Product Information

Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	RM153
Preservative	0.05% Sodium Azide
Isotype	IgG Kappa
Conjugate	mFluor Violet 450 SE
Purity	Protein A purified
Buffer	50mM Sodium Borate

Product Description

Host	Rabbit
Gene ID	944
Gene Symbol	TNFSF8
Species	Mouse
Specificity/Sensitivity	This antibody is specific for murine CD30 Ligand/TNFSF8, a 40-kDa type II transmembrane protein member of the TNF superfamily. Interaction of CD30 with CD30 Ligand/TNFSF8 plays an important role in T cell activation/costimulation.
Immunogen	This antibody was raised by immunising Sprague Dawley rats with the CHO cells stably transfected with murine CD30 Ligand/TNFSF8 and fusing the splenocytes with P3U1 myeloma cells.

Product Application Details

Applications	Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Block/Neutralize, Functional Assay
Recommended Dilutions	Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Block/Neutralize, Functional Assay
Application Notes	Optimal dilution of this antibody should be experimentally determined.

Images

CD30 Ligand/TNFSF8 Antibody (RM153) [mFluor Violet 450 SE] [NBP3-11991MFV450] - Vial of mFluor Violet 450 conjugated antibody. mFluor Violet 450 is optimally excited at 406 nm by the Violet laser (405 nm) and has an emission maximum of 445 nm.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-11991MFV450

210-TA-005	TNF-alpha [Unconjugated]
1028-CL-050	CD30 Ligand/TNFSF8 [Unconjugated]
M6000B-1	IL-6 [HRP]
6507-IL-010/CF	IL-4 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-11991MFV450

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

