

Product Datasheet

IL-9 Antibody (RM9A4) [mFluor Violet 450 SE] NBP3-12023MFV450

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-12023MFV450

Updated 7/11/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-12023MFV450



NBP3-12023MFV450

IL-9 Antibody (RM9A4) [mFluor Violet 450 SE]

| Product Information | |
|---------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | RM9A4 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG Kappa |
| Conjugate | mFluor Violet 450 SE |
| Purity | Protein A purified |
| Buffer | 50mM Sodium Borate |

| Product Description | |
|-------------------------|---|
| Description | mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet. |
| Host | Rabbit |
| Gene ID | 3578 |
| Gene Symbol | IL9 |
| Species | Mouse |
| Specificity/Sensitivity | This antibody is specific for murine IL-9. |
| Immunogen | This antibody was raised by immunising rats with recombinant full-length murine IL-9. |

| Product Application Details | |
|-----------------------------|--|
| Applications | Flow Cytometry |
| Recommended Dilutions | Flow Cytometry |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |

Images

IL-9 Antibody (RM9A4) [mFluor Violet 450 SE] [NBP3-12023MFV450] - Vial of mFluor Violet 450 conjugated antibody. mFluor Violet 450 is optimally excited at 406 nm by the Violet laser (405 nm) and has an emission maximum of 445 nm.



mFluor™ Violet 450

| | |
|---------------------|-------------------|
| LASER (nm) | FILTER |
| Violet (405) | 450/45 |
| EXCITATION MAX (nm) | EMISSION MAX (nm) |
| 406 | 445 |



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to NBP3-12023MFV450

| | |
|----------------|------------------------------|
| NBP2-35257-2ug | Recombinant Rat IL-9 Protein |
| 210-TA-005 | TNF-alpha [Unconjugated] |
| 409-ML-010/CF | IL-9 [Unconjugated] |
| D6050 | IL-6 [HRP] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-12023MFV450

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

