

# Product Datasheet

## MS4A4A Antibody (3F2) [DyLight 488] NBP3-12045G

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-12045G](http://www.novusbio.com/NBP3-12045G)

Updated 7/11/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-12045G](http://www.novusbio.com/reviews/destination/NBP3-12045G)



**NBP3-12045G**

MS4A4A Antibody (3F2) [DyLight 488]

| Product Information         |  |
|-----------------------------|--|
| Unit Size                   | 0.1 ml   |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.  |
| Storage                     | Store at 4C in the dark.   |
| Clonality                   | Monoclonal   |
| Clone                       | 3F2  |
| Preservative                | 0.05% Sodium Azide   |
| Isotype                     | IgG Kappa  |
| Conjugate                   | DyLight 488  |
| Purity                      | Protein A purified   |
| Buffer                      | 50mM Sodium Borate   |
| Product Description         |  |
| Host                        | Rabbit   |
| Gene ID                     | 51338  |
| Gene Symbol                 | MS4A4A   |
| Species                     | Human  |
| Specificity/Sensitivity     | This antibody is specific for human extracellular MS4A4A. MS4A (membrane-spanning 4-domain family, subfamily A) is a large family of proteins that includes at least 26 members in mouse and humans. Flanked by amino- and carboxyl-cytoplasmic regions, MS4A family members contain four highly conserved transmembrane domains connected by short extracellular and cytoplasmic loops. MS4A4A is likely to be a component of oligomeric cell surface complexes involved in signal transduction in diverse cell lineages. |
| Immunogen                   | This antibody was generated by immunizing C57BL/6 mice against mouse cell line 300-19 expressing transfected human membrane-spanning 4-domains subfamily A member 4A (MS4A4A)-GFP.   |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.  |
| Product Application Details |  |
| Applications                | Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation   |
| Recommended Dilutions       | Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation   |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.   |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-12045G**

---

|               |                                    |
|---------------|------------------------------------|
| NBP1-86542PEP | MS4A4A Recombinant Protein Antigen |
| 7268-CT-100   | CTLA-4 [Unconjugated]              |
| NBL1-13318    | MS4A4A Overexpression Lysate       |
| NBP2-22203    | ERK1 Antibody (1E5)                |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-12045G](http://www.novusbio.com/reviews/submit/NBP3-12045G)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

