

Product Datasheet

UCN2 Antibody NBP3-12450

Unit Size: 100 ug

Store at 4C short term. Store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-12450

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-12450



NBP3-12450

UCN2 Antibody

Product Information

Unit Size	100 ug
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Purity	Affinity purified
Buffer	Contains Tris, HCl/glycine buffer (pH 7.4-7.8), 30% glycerol, 0.5% BSA, along with cryo-protective agents, and HEPES (exact storage and stabilization buffer information is proprietary).

Product Description


Description	Novus Biologicals Rabbit UCN2 Antibody (NBP3-12450) is a polyclonal antibody validated for use in IHC, WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	90226
Gene Symbol	UCN2
Species	Human, Mouse
Immunogen	Synthetic peptide taken within amino acid region 70-115 of human Urocortin 2 protein.

Product Application Details


Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500, ELISA 1:10000, Immunohistochemistry, Immunohistochemistry-Paraffin

Images


Immunohistochemistry-Paraffin: UCN2 Antibody [NBP3-12450] - Mouse liver. 1:1000 dilution in IHC blocking buffer. DAB (brown) staining in Hematoxylin QS (blue) counterstain. 40X magnification.



Immunohistochemistry-Paraffin: UCN2 Antibody [NBP3-12450] - Mouse colon. 1:1000 dilution in IHC blocking buffer. DAB (brown) staining in Hematoxylin QS (blue) counterstain. 5x magnification.




Immunohistochemistry-Paraffin: UCN2 Antibody [NBP3-12450] - Mouse colon. 1:1000 dilution in IHC blocking buffer. DAB (brown) staining in Hematoxylin QS (blue) counterstain. 40 magnification.



Immunohistochemistry-Paraffin: UCN2 Antibody [NBP3-12450] - Mouse liver. 1:1000 dilution in IHC blocking buffer. DAB (brown) staining in Hematoxylin QS (blue) counterstain. 5X magnification.



Southern Blot: UCN2 Antibody [NBP3-12450] - 1:500 antibody dilution in DiluObuffer.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-12450

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-57622PEP	UCN2 Recombinant Protein Antigen
9134-TN-050	C1QL1 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-12450

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

