

Product Datasheet

EBOV GP Antibody - Azide and BSA Free NBP3-12747-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-12747

Updated 3/3/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-12747



NBP3-12747-100ul

EBOV GP Antibody - Azide and BSA Free


Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Antigen and protein A Affinity-purified
Buffer	0.2 um filtered solution in PBS

Product Description	
Description	This antibody can be stored at 2-8C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20C to -80C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Gene ID	3160774
Gene Symbol	GP
Species	Ebola Virus
Reactivity Notes	Ebola virus EBOV (subtype Zaire, strain Mayinga 1976) Has cross-reactivity in WB with: Ebola virus EBOV (subtype Zaire, strain H.sapiens-wt/GIN/2014/Kissidougou-C15) Glycoprotein GP Protein (His Tag) Ebola virus EBOV (subtype Bundibugyo, strain Uganda 2007) Glycoprotein GP Protein Ebola virus EBOV (Subtype Sudan, strain Gulu) Glycoprotein GP Protein (aa:Met1-Asn637, His Tag)
Specificity/Sensitivity	Ebola virus EBOV (subtype Zaire, strain Mayinga 1976) Glycoprotein GP
Immunogen	Produced in rabbits immunized with purified, recombinant EBOV GP (Accession#: AAC54887.1; Met1-Gln650)

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot 1:500-1:2000

Images

Western Blot: EBOV GP Antibody [NBP3-12747] - Anti-Ebola virus EBOV (subtype Zaire, strain Mayinga 1976) Glycoprotein / GP rabbit polyclonal antibody at 1:5000 dilution Lane A:Ebola virus EBOV (subtype Zaire, strain H.sapiens-wt/GIN/2014/Kissidougou-C15) Glycoprotein / GP at 10ng Lane B:GP (B-EBOV-2007) at 10ng Lane C:Ebola virus EBOV (subtype Zaire, strain Mayinga 1976) Glycoprotein / GP at 10ng Lane D:Ebola virus EBOV (Subtype Sudan, strain Gulu) Glycoprotein / GP Protein (aa:Met1-Asn637) at 10ng Secondary HRP Conjugated Goat-anti-Rabbit IgG (H+L) at 1/1000 dilution Performed under reducing conditions.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950

Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400

Toll Free: 855.668.8722

Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15

Fax: (44) (0) 1235 533420

info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com

Technical Support: nb-technical@bio-techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Products Related to NBP3-12747-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
9016-EB-100	EBOV GP [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-12747

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

