

Product Datasheet

EBOV VP40 protein Antibody - Azide and BSA Free NBP3-12760-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-12760

Updated 3/3/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-12760



NBP3-12760-100ul

EBOV VP40 protein Antibody - Azide and BSA Free


Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Antigen and protein A Affinity-purified
Buffer	0.2 um filtered solution in PBS

Product Description	
Description	This antibody can be stored at 2-8C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20C to -80C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Species	Ebola Virus
Reactivity Notes	Ebola virus EBOV (subtype Zaire, strain H.sapiens-wt GIN 2014 Kissidougou-C15)
Specificity/Sensitivity	Ebola virus EBOV (subtype Zaire, strain H.sapiens-wt GIN 2014 Kissidougou-C15) Matrix protein VP40
Immunogen	Produced in rabbits immunized with purified, recombinant EBOV VP40 protein (Accession#: AHX24648.1; Asn31-Lys326)

Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 1:1000-1:5000, ELISA 1:5000-1:10000

Images

Western Blot: EBOV VP40 protein Antibody [NBP3-12760] - nti-Ebola virus EBOV (subtype Zaire, strain H.sapiens-wt/GIN/2014/Kissidougou-C15) Matrix protein VP40 rabbit polyclonal antibody at 1:5000 dilution
 Lane A: Ebola virus EBOV (subtype Sudan, strain Gulu) VP40 / Matrix protein VP40 Protein (10ng) Lane B: Ebola virus EBOV (subtype Sudan, strain Gulu) VP40 / Matrix protein VP40 Protein (10ng) Lane C: Ebola virus EBOV (subtype Bundibugyo, strain Uganda 2007) VP40 (10ng) Lane D: Ebola virus EBOV (subtype Zaire, strain H.sapiens-wt/GIN/2014/Kissidougou-C15) VP40 (10ng) Lane E: Ebola virus EBOV (subtype Zaire, strain H.sapiens-wt/GIN/2014/Kissidougou-C15) Matrix protein VP40 Protein(10ng) Lane F: Human CLIC1(10ng) Secondary Dylight 800-labeled Antibody To Rabbit IgG (H+L) at 1:5000 dilution. Developed using the Odyssey technique. Performed under reducing conditions.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-12760-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-12760

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

