

Product Datasheet

Influenza A H3N2 Hemagglutinin Antibody - (A/Taiwan/3446/02) NBP3-12877

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-12877

Updated 5/1/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-12877



NBP3-12877

Influenza A H3N2 Hemagglutinin Antibody - (A/Taiwan/3446/02)

Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.025% Proclin 300
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS (pH 7), 20% Glycerol


Product Description	
Description	Centrifuge briefly prior to opening.
Host	Rabbit
Species	Influenza A Virus H3N2
Specificity/Sensitivity	Influenza A H3N2 Hemagglutinin (A/Taiwan/3446/02)
Immunogen	Recombinant protein encompassing a sequence within the N-terminus region of Influenza A H3N2 Hemagglutinin (Hemagglutinin) (A/TW/3446/02(H3N2). The exact sequence is proprietary.

Product Application Details	
Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin



Images


Western Blot: Influenza A H3N2 Hemagglutinin Antibody - (A/Taiwan/3446/02) [NBP3-12877] - Recombinant H3N2 HA protein were separated by 10% SDS-PAGE, and the membrane was blotted with Influenza A H3N2 Hemagglutinin (Hemagglutinin) antibody (NBP3-12877) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.



Immunohistochemistry: Influenza A H3N2 Hemagglutinin Antibody - (A/Taiwan/3446/02) [NBP3-12877] - Influenza A H3N2 Hemagglutinin antibody detects Influenza A Virus H3N2 Hemagglutinin (HA) protein on mouse nasal tissue (infect with Influenza Virus A) by immunohistochemical analysis. Sample: Paraffin-embedded mouse nasal tissue (infect with Influenza Virus A). Influenza A Virus H3N2 Hemagglutinin (HA) antibody (NBP3-12877) dilution: 1:50. Antigen Retrieval: Trilogy(TM) (EDTA based, pH 8.0) buffer, 15min



Western Blot: Influenza A H3N2 Hemagglutinin Antibody - (A/Taiwan/3446/02) [NBP3-12877] - Influenza A Virus H3N2 Hemagglutinin (HA) antibody detects HA (H3N2) protein by western blot analysis. A. MDCK cell lysates infected with 4 HAU/ml of A/duck/Alberta/35/1976-like, H1N1 B. MDCK cell lysates infected with 4 HAU/ml of A/duck/Ukraine/1963-like, H3N8 C. MDCK cell lysates infected with 4 HAU/ml of A/mallard/Miyagi/53/1976-like, H5N3 D. MDCK cell lysates infected with 4 HAU/ml of A/duck/Taiwan/33/1993-like, H7N7 E. MDCK cell lysates infected with 4 HAU/ml of A/turkey/Wisconsin/1/1966-like, H9N2 F. MDCK cell lysates infected with 4 HAU/ml of A/PR/8/34, H1N1 G. MDCK cell lysates infected with 4 HAU/ml of A/Taiwan/126/2009, pdmH1N1 H. MDCK cell lysates infected with 4 HAU/ml of A/Taiwan/3446/2002, H3N2 12% SDS-PAGE Influenza A Virus H3N2 Hemagglutinin (HA) antibody (NBP3-12877) dilution: 1:1000 The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-12877

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-99037-100ug	Recombinant Influenza A H3N2 Hemagglutinin (A/canine/Guangdong/1/2006) His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-12877

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

