

Product Datasheet

MBL-2 Antibody NBP3-13134

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-13134

Updated 5/1/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-13134




NBP3-13134**MBL-2 Antibody**

| Product Information | |
|------------------------------------|---|
| Unit Size | 100 ul |
| Concentration | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| Clonality | Polyclonal |
| Preservative | 0.025% Proclin 300 |
| Isotype | IgG |
| Purity | Antigen Affinity-purified |
| Buffer | PBS (pH 7), 20% Glycerol, 1% BSA |
| Product Description | |
| Description | Centrifuge briefly prior to opening. |
| Host | Rabbit |
| Gene ID | 17195 |
| Gene Symbol | Mbl2 |
| Species | Human, Mouse |
| Reactivity Notes | Rhesus Monkey (89%) |
| Immunogen | Recombinant protein encompassing a sequence within the center region of human MBL-2. The exact sequence is proprietary. |
| Product Application Details | |
| Applications | Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry |
| Recommended Dilutions | Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin |




Images


Western Blot: MBL-2 Antibody [NBP3-13134] - Human plasma (30 ug) was separated by 12% SDS-PAGE, and the membrane was blotted with MBL-2 antibody (NBP3-13134) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.



Immunohistochemistry-Paraffin: MBL-2 Antibody [NBP3-13134] - MBL-C antibody detects MBL-C protein at cytoplasm by immunohistochemical analysis. Sample: Paraffin-embedded mouse liver. MBL-C stained by MBL-C antibody (NBP3-13134) diluted at 1:500. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



Western Blot: MBL-2 Antibody [NBP3-13134] - Various tissue extracts (50 ug) were separated by 12% SDS-PAGE, and the membrane was blotted with MBL-C antibody (NBP3-13134) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-13134

| | |
|-----------------|---|
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |
| NBP2-61994-20ug | Recombinant Human MBL-2 His Protein |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-13134

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

