

Product Datasheet

Influenza A H1N1 Hemagglutinin Antibody (GT223) - (A/WSN/1933) NBP3-13497

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-13497

Updated 6/14/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-13497



NBP3-13497

Influenza A H1N1 Hemagglutinin Antibody (GT223) - (A/WSN/1933)

Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	GT223
Preservative	No Preservative
Isotype	IgG1
Purity	Protein G purified
Buffer	PBS


Product Description	
Description	Centrifuge briefly prior to opening.
Host	Mouse
Species	Influenza A Virus H1N1
Specificity/Sensitivity	Influenza A H1N1 Hemagglutinin. (A/WSN/1933)
Immunogen	Recombinant protein encompassing a sequence within the C-terminus region of Influenza A Virus H1N1 Neuraminidase (Hemagglutinin) (A/WSN/1933(H1N1)). The exact sequence is proprietary.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot, Immunocytochemistry/ Immunofluorescence




Images

Western Blot: Influenza A H1N1 Hemagglutinin Antibody (GT223) - (A/WSN/1933) [NBP3-13497] - Influenza A H1N1 Hemagglutinin (HA) antibody [GT223] detects Influenza A Virus H1N1 Hemagglutinin (HA) protein by Western blot analysis. A. 5 ug DF1 whole cell lysate/extract B. 5 ug whole cell lysate/extract of Influenza A (WSN) infected DF1 cells 7.5 % SDS-PAGE Influenza A Virus H1N1 Hemagglutinin (HA) antibody [GT223] (NBP3-13497) dilution: 1:1000



Immunocytochemistry/Immunofluorescence: Influenza A H1N1 Hemagglutinin Antibody (GT223) - (A/WSN/1933) [NBP3-13497] - Influenza A H1N1 Hemagglutinin (HA) antibody [GT223] detects Influenza A Virus H1N1 Hemagglutinin (HA) protein by immunofluorescent analysis. Sample: DF1 cells infected with Influenza A virus (WSN) . Green: Influenza A Virus H1N1 Hemagglutinin (HA) protein stained by Influenza A Virus H1N1 Hemagglutinin (HA) antibody [GT223] (NBP3-13497) diluted at 1:500. Blue: Hoechst 33342 staining.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-13497

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NBP1-99041-50ug	Recombinant Influenza A H1N1 Hemagglutinin (A/New York/3571/2009) His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-13497

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

