

Product Datasheet

T Plastin Antibody (GT3310) NBP3-13518

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-13518

Updated 3/7/2022 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-13518



NBP3-13518

T Plastin Antibody (GT3310)


Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	GT3310
Preservative	No Preservative
Isotype	IgG2b
Purity	Protein A purified
Buffer	PBS

Product Description	
Host	Mouse
Gene ID	5358
Gene Symbol	PLS3
Species	Human, Mouse, Rat
Reactivity Notes	Zebrafish (89%)
Immunogen	Recombinant protein encompassing a sequence within the center region of human T Plastin. The exact sequence is proprietary.
Notes	Centrifuge briefly prior to opening.


Product Application Details	
Applications	Western Blot, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin

Images


Western Blot: T Plastin Antibody (GT3310) [NBP3-13518] - Non-transfected (-) and transfected (+) 293T whole cell extracts (30 ug) were separated by 7.5% SDS-PAGE, and the membrane was blotted with T Plastin antibody [GT3310] (NBP3-13518) diluted at 1:1000.




Immunocytochemistry/Immunofluorescence: T Plastin Antibody (GT3310) [NBP3-13518] - T Plastin antibody [GT3310] detects T Plastin protein at cytoplasm by immunofluorescent analysis. Sample: HeLa cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: T Plastin protein stained by T Plastin antibody [GT3310] (NBP3-13518) diluted at 1:200. Blue: Hoechst 33342 staining. Scale bar = 10 μ m.




Immunohistochemistry-Paraffin: T Plastin Antibody (GT3310) [NBP3-13518] - T Plastin antibody [GT3310] detects T Plastin protein at cytoplasm in rat liver by immunohistochemical analysis. Sample: Paraffin-embedded rat liver. T Plastin antibody [GT3310] (NBP3-13518) diluted at 1:250. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



Western Blot: T Plastin Antibody (GT3310) [NBP3-13518] - T Plastin antibody detects T Plastin protein by western blot analysis. Various whole cell extracts (30 μ g) were separated by 7.5% SDS-PAGE, and the membrane was blotted with T Plastin antibody (NBP3-13518) diluted at a dilution of 1:1000.



Immunohistochemistry-Paraffin: T Plastin Antibody (GT3310) [NBP3-13518] - T Plastin antibody [GT3310] detects T Plastin protein at cytoplasm in mouse cervix by immunohistochemical analysis. Sample: Paraffin-embedded mouse cervix. T Plastin antibody [GT3310] (NBP3-13518) diluted at 1:250. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to NBP3-13518

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP2-27231	Mouse IgG2b Isotype Control (MPC-11)
H00005358-P01-10ug	Recombinant Human T Plastin GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-13518

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

