

Product Datasheet

EEN Antibody (GT3065) NBP3-13631

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-13631

Updated 3/8/2022 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-13631



NBP3-13631

EEN Antibody (GT3065)

Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	GT3065
Preservative	No Preservative
Isotype	IgG2a
Purity	Protein A purified
Buffer	PBS


Product Description	
Host	Mouse
Gene ID	6455
Gene Symbol	SH3GL1
Species	Human, Mouse, Rat
Reactivity Notes	Human (85%), Bovine (83%)
Immunogen	Recombinant protein encompassing a sequence within the center region of mouse EEN. The exact sequence is proprietary.
Notes	Centrifuge briefly prior to opening.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin




Images


Western Blot: EEN Antibody (GT3065) [NBP3-13631] - Various tissue extracts (50 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with EEN antibody [GT3065] (NBP3-13631) diluted at 1:1000. The HRP-conjugated anti-mouse IgG antibody (NBP2-19382) was used to detect the primary antibody.



Immunocytochemistry/Immunofluorescence: EEN Antibody (GT3065) [NBP3-13631] - EEN antibody [GT3065] detects EEN protein at cytoplasm by immunofluorescent analysis. Sample: HeLa cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: EEN stained by EEN antibody [GT3065] (NBP3-13631) diluted at 1:500. Blue: Hoechst 33342 staining.



Immunohistochemistry-Paraffin: EEN Antibody (GT3065) [NBP3-13631] - EEN antibody [GT3065] detects EEN protein at cytoplasm by immunohistochemical analysis. Sample: Paraffin-embedded mouse hippocampus. EEN stained by EEN antibody [GT3065] (NBP3-13631) diluted at 1:200. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to NBP3-13631

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96778	Mouse IgG2a Isotype Control (M2A)
H00006455-Q01-10ug	Recombinant Human EEN GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-13631

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

