

Product Datasheet

Serum Response Factor SRF Antibody (PCRP-SRF-1F1) NBP3-13784-100ug

Unit Size: 100 ug

Store at 4C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-13784

Updated 7/16/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-13784



NBP3-13784-100ug

Serum Response Factor SRF Antibody (PCRP-SRF-1F1)

Product Information	
Unit Size	100 ug
Concentration	0.2 mg/ml
Storage	Store at 4C.
Clonality	Monoclonal
Clone	PCRP-SRF-1F1
Preservative	0.05% Sodium Azide
Isotype	IgG1
Purity	Protein A or G purified
Buffer	10 mM PBS with 0.05% BSA


Product Description	
Description	<p>200ug/ml of antibody purified from Bioreactor Concentrate by Protein A or G. Prepared in 10 mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0 mg/ml. (NBP3-14187)</p> <p>Antibody with azide - store at 2 to 8C. Antibody without azide - store at -20 to -80 C. Non-hazardous.</p> <p>Theoretical molecular weight: 40-67kDa</p>
Host	Mouse
Gene ID	6722
Gene Symbol	SRF
Species	Human
Immunogen	Recombinant full-length human Serum Response Factor SRF protein (Uniprot: P11831)

Product Application Details	
Applications	Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation, Protein Array
Recommended Dilutions	Western Blot 1-2 ug/ml, Flow Cytometry 1-2 ug/million cells, ELISA, Immunocytochemistry/ Immunofluorescence 1-2 ug/ml, Immunoprecipitation 1-2 ug/100-500 ug protein, Protein Array
Application Notes	<p>ELISA: For coating, order antibody without BSA</p> <p>Immunoprecipitation: 1-2ug per 100-500ug of total protein (1ml of cell lysate)</p>




Images


Immunocytochemistry/Immunofluorescence: Serum Response Factor SRF Antibody (PCRP-SRF-1F1) [NBP3-13784] - Immunofluorescence Analysis of PFA-fixed HeLa cells stained using Serum Response Factor SRF Antibody (PCRP-SRF-1F1) followed by goat anti-mouse IgG-CF488 (green). CF640R phalloidin (red).



Flow Cytometry: Serum Response Factor SRF Antibody (PCRP-SRF-1F1) [NBP3-13784] - Flow cytometric analysis of PFA-fixed HeLa cells. Serum Response Factor SRF antibody (PCRP-SRF-1F1) followed by goat anti-mouse IgG-CF488 (blue); isotype control (red).



Protein Array: Serum Response Factor SRF Antibody (PCRP-SRF-1F1) [NBP3-13784] - Analysis of Protein Array containing more than 19,000 full-length human proteins using Serum Response Factor SRF Antibody (PCRP-SRF-1F1).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-13784-100ug

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
H00006722-Q01-10ug	Recombinant Human Serum Response Factor SRF GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-13784

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

