

Product Datasheet

SLC45A3/Prostein Antibody (SLC45A3/7176R) NBP3-13804-100ug

Unit Size: 100 ug

Store at 4C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-13804

Updated 7/16/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-13804



NBP3-13804-100ug

SLC45A3/Prostein Antibody (SLC45A3/7176R)

| Product Information | |
|-------------------------|--------------------------|
| Unit Size | 100 ug |
| Concentration | 0.2 mg/ml |
| Storage | Store at 4C. |
| Clonality | Monoclonal |
| Clone | SLC45A3/7176R |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG Kappa |
| Purity | Protein A or G purified |
| Buffer | 10 mM PBS with 0.05% BSA |
| Target Molecular Weight | 59 kDa |


| Product Description | |
|---------------------|---|
| Description | 200ug/ml of antibody purified from Bioreactor Concentrate by Protein A or G. Prepared in 10 mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0 mg/ml. (NBP3-14049) Antibody with azide - store at 2 to 8C. Antibody without azide - store at -20 to -80 C. |
| Host | Rabbit |
| Gene ID | 85414 |
| Gene Symbol | SLC45A3 |
| Species | Human |
| Immunogen | Synthesized peptides to the N terminus of human SLC45A3/Prostein protein (Uniprot: Q96JT2) |

| Product Application Details | |
|-----------------------------|--|
| Applications | Immunohistochemistry-Paraffin, Protein Array |
| Recommended Dilutions | Immunohistochemistry-Paraffin 1-2 ug/ml, Protein Array |
| Application Notes | Immunohistochemistry (Formalin-fixed): 1-2ug/ml for 30 minutes at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95C followed by cooling at RT for 20 minutes |



Images

Immunohistochemistry-Paraffin: SLC45A3/Prostein Antibody (SLC45A3/7176R) [NBP3-13804] - Formalin-fixed, paraffin-embedded human prostate stained with SLC45A3/Prostein antibody (SLC45A3/7176R). Inset: PBS instead of primary antibody; secondary only negative control.



Protein Array: SLC45A3/Prostein Antibody (SLC45A3/7176R) [NBP3-13804] - Analysis of Protein Array containing more than 19,000 full-length human proteins using SLC45A3/Prostein Antibody (SLC45A3/7176R).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-13804-100ug

| | |
|---------------|---|
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP1-89629PEP | SLC45A3/Prostein Recombinant Protein Antigen |
| DKK300 | Kallikrein 3/PSA [HRP] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-13804

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

