

# Product Datasheet

## GLIS3 Antibody (PCRP-GLIS3-1B11) [DyLight 488] NBP3-14052G

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-14052G](http://www.novusbio.com/NBP3-14052G)

Updated 10/26/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-14052G](http://www.novusbio.com/reviews/destination/NBP3-14052G)



**NBP3-14052G**

GLIS3 Antibody (PCRP-GLIS3-1B11) [DyLight 488]

| Product Information         |                                                                                             |
|-----------------------------|---------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                      |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                     | Store at 4C in the dark.                                                                    |
| Clonality                   | Monoclonal                                                                                  |
| Clone                       | PCRP-GLIS3-1B11                                                                             |
| Preservative                | 0.05% Sodium Azide                                                                          |
| Isotype                     | IgG1 Kappa                                                                                  |
| Conjugate                   | DyLight 488                                                                                 |
| Purity                      | Protein A or G purified                                                                     |
| Buffer                      | 50mM Sodium Borate                                                                          |
| Product Description         |                                                                                             |
| Host                        | Mouse                                                                                       |
| Gene ID                     | 169792                                                                                      |
| Gene Symbol                 | GLIS3                                                                                       |
| Species                     | Human                                                                                       |
| Immunogen                   | Recombinant full-length human GLIS3 protein (Uniprot: Q8NEA6)                               |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.           |
| Product Application Details |                                                                                             |
| Applications                | ELISA, Flow Cytometry, Immunoprecipitation, Protein Array                                   |
| Recommended Dilutions       | Flow Cytometry, ELISA, Immunoprecipitation, Protein Array                                   |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                      |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NBP3-14052G**

---

|                   |                                                             |
|-------------------|-------------------------------------------------------------|
| NBP1-43319G-0.5ml | Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1) [DyLight 488] |
| NBP2-33787PEP     | GLIS3 Recombinant Protein Antigen                           |
| 291-G1-200        | IGF-I/IGF-1 [Unconjugated]                                  |
| NBP2-05335        | GLIS3 Overexpression Lysate                                 |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-14052G](http://www.novusbio.com/reviews/submit/NBP3-14052G)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

