

# Product Datasheet

## IRF3 Antibody (PCRP-IRF3-3B2) [DyLight 550] NBP3-14106R

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-14106R](http://www.novusbio.com/NBP3-14106R)

Updated 10/26/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-14106R](http://www.novusbio.com/reviews/destination/NBP3-14106R)



**NBP3-14106R**

IRF3 Antibody (PCRP-IRF3-3B2) [DyLight 550]

| Product Information |   |
|---------------------|---|
| Unit Size           | 0.1 ml  |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.  |
| Clonality           | Monoclonal  |
| Clone               | PCRP-IRF3-3B2   |
| Preservative        | 0.05% Sodium Azide  |
| Isotype             | IgG1 Kappa  |
| Conjugate           | DyLight 550   |
| Purity              | Protein A or G purified   |
| Buffer              | 50mM Sodium Borate  |

| Product Description |   |
|---------------------|---|
| Host                | Mouse   |
| Gene ID             | 3661  |
| Gene Symbol         | IRF3  |
| Species             | Human   |
| Reactivity Notes    | 0   |
| Immunogen           | Recombinant full-length human IRF3 protein (Uniprot: Q14653)                      |
| Notes               | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |

| Product Application Details |  |
|-----------------------------|--|
| Applications                | Western Blot, Immunocytochemistry/Immunofluorescence, Protein Array    |
| Recommended Dilutions       | Western Blot, Immunocytochemistry/Immunofluorescence, Protein Array    |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NBP3-14106R**

---

|             |   |
|-------------|---|
| NBP1-43319R | Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1) [DyLight 550] |
| NBC1-18469  | Recombinant Human IRF3 Protein                              |
| 210-TA-005  | TNF-alpha [Unconjugated]                                    |
| NBL1-12032  | IRF3 Overexpression Lysate                                  |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-14106R](http://www.novusbio.com/reviews/submit/NBP3-14106R)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

