

# Product Datasheet

## **FABP3/H-FABP Antibody (FABP3/3430) [DyLight 550] NBP3-14215R**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-14215R](http://www.novusbio.com/NBP3-14215R)

Updated 10/26/2023 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-14215R](http://www.novusbio.com/reviews/destination/NBP3-14215R)



**NBP3-14215R**

FABP3/H-FABP Antibody (FABP3/3430) [DyLight 550]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | FABP3/3430                                                                                  |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG2b Kappa                                                                                 |
| Conjugate           | DyLight 550                                                                                 |
| Purity              | Protein A or G purified                                                                     |
| Buffer              | 50mM Sodium Borate                                                                          |

| Product Description |                                                                                                                    |
|---------------------|--------------------------------------------------------------------------------------------------------------------|
| Host                | Mouse                                                                                                              |
| Gene ID             | 2170                                                                                                               |
| Gene Symbol         | FABP3                                                                                                              |
| Species             | Human                                                                                                              |
| Reactivity Notes    | 0                                                                                                                  |
| Immunogen           | Human recombinant FABP3/H-FABP protein fragment (around aa1-127) (exact sequence is proprietary) (Uniprot: P05413) |
| Notes               | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                  |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Immunohistochemistry-Paraffin, Protein Array                           |
| Recommended Dilutions       | Immunohistochemistry-Paraffin, Protein Array                           |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NBP3-14215R**

---

|                |                                                                    |
|----------------|--------------------------------------------------------------------|
| NBP1-43317R    | Mouse IgG2b Kappa Light Chain Isotype Control (MG2b) [DyLight 550] |
| NBP3-07107-5ug | Recombinant Human FABP3/H-FABP Protein                             |
| DY1857-05      | Lipocalin-2/NGAL [Biotin]                                          |
| DY1678         | FABP3/H-FABP [Biotin]                                              |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-14215R](http://www.novusbio.com/reviews/submit/NBP3-14215R)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

