

# Product Datasheet

## TADA1L Antibody (PCRP-TADA1-1C9) [mFluor Violet 500 SE] NBP3-14264MFV500

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-14264MFV500](http://www.novusbio.com/NBP3-14264MFV500)

Updated 10/26/2023 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-14264MFV500](http://www.novusbio.com/reviews/destination/NBP3-14264MFV500)



**NBP3-14264MFV500**

TADA1L Antibody (PCRP-TADA1-1C9) [mFluor Violet 500 SE]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	PCRP-TADA1-1C9
Preservative	0.05% Sodium Azide
Isotype	IgG2b
Conjugate	mFluor Violet 500 SE
Purity	Protein A or G purified
Buffer	50mM Sodium Borate

Product Description	
Host	Mouse
Gene ID	117143
Gene Symbol	TADA1
Species	Human
Reactivity Notes	Predicted to react with Rabbit, Goat, Equine, Guinea Pig, Canine, Porcine, Rhesus Monkey, Orangutan and Elephant (based on sequence homology).
Immunogen	Recombinant full-length human TADA1Lprotein (Uniprot: Q96BN2)
Notes	mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.

Product Application Details	
Applications	Western Blot, Flow Cytometry, Immunoprecipitation, Protein Array
Recommended Dilutions	Western Blot, Flow Cytometry, Immunoprecipitation, Protein Array
Application Notes	Optimal dilution of this antibody should be experimentally determined.

**Images**

TADA1L Antibody (PCRP-TADA1-1C9) [mFluor Violet 500 SE] - Vial of mFluor Violet 500 conjugated antibody. mFluor Violet 500 is optimally excited at 410 nm by the Violet laser (405 nm) and has an emission maximum of 501 nm.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NBP3-14264MFV500**

---

H00117143-P01-2ug	Recombinant Human TADA1L GST (N-Term) Protein
2695-SE-010	Airway Trypsin-like Protease/HAT/TMPRSS11D
AF1180	DPPIV/CD26 Antibody [Unconjugated]
NB110-41404	Glutamine Synthetase Antibody

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-14264MFV500](http://www.novusbio.com/reviews/submit/NBP3-14264MFV500)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

