

# Product Datasheet

## Human Coronavirus Spike Protein Antibody (15) [DyLight 594]

**NBP3-14656DL594**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-14656DL594](http://www.novusbio.com/NBP3-14656DL594)

Updated 4/19/2022 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-14656DL594](http://www.novusbio.com/reviews/destination/NBP3-14656DL594)



**NBP3-14656DL594**

Human Coronavirus Spike Protein Antibody (15) [DyLight 594]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 15                                                                                          |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG1                                                                                        |
| Conjugate           | DyLight 594                                                                                 |
| Purity              | Protein A purified                                                                          |
| Buffer              | 50mM Sodium Borate                                                                          |

| Product Description     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|-------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host                    | Mouse                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Gene ID                 | 3200426                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Gene Symbol             | S                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Species                 | HCoV-OC43                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Specificity/Sensitivity | <p>Recognizes the Spike of HCoV-OC43. Doesn't recognize the Spike of SARS-CoV-2, SARS-CoV, MERS-CoV, HCoV-NL63, HCoV-229E, HCoV-HKU1 (isolate N5).</p> <p>Has cross-reactivity in ELISA and WB with:<br/>HCoV-OC43 Spike S1+S2 ECD-His</p> <p>No cross-reactivity in ELISA and WB with:<br/>SARS-CoV-2 Spike S1-His<br/>SARS-CoV-2 Spike S1 NTD-Fc &amp; AVI<br/>SARS-CoV-2 Spike S1 NTD-His &amp; AVI<br/>SARS-CoV-2 Spike RBD-His<br/>SARS-CoV-2 Spike S1+S2 ECD-His<br/>SARS-CoV-2 Spike S2 ECD-His<br/>SARS-CoV Spike S1+S2 ECD-His<br/>MERS-CoV Spike S1+S2 ECD-His<br/>HCoV-229E Spike S1+S2 ECD-His<br/>HCoV-HKU1 (isolate N5) Spike S1+S2 ECD-His<br/>HCoV-NL63 Spike S1+S2 ECD-His</p> |
| Immunogen               | This product is a recombinant monoclonal antibody expressed from HEK293 cells.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Notes                   | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Western Blot, ELISA                                                    |
| Recommended Dilutions       | Western Blot, ELISA                                                    |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NBP3-14656DL594**

---

|                 |                                                |
|-----------------|------------------------------------------------|
| NBP1-97005DL594 | Mouse IgG1 Isotype Control (MG1) [DyLight 594] |
|-----------------|------------------------------------------------|

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-14656DL594](http://www.novusbio.com/reviews/submit/NBP3-14656DL594)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

