

Product Datasheet

CD27 Ligand/TNFSF7/CD70 Antibody (BLR176J) NBP3-14748

Unit Size: 100 ul

Store at 2 - 8 C / 1 year from date of receipt

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-14748

Updated 9/12/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-14748



NBP3-14748

CD27 Ligand/TNFSF7/CD70 Antibody (BLR176J)

Product Information

Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 2 - 8 C / 1 year from date of receipt
Clonality	Monoclonal
Clone	BLR176J
Preservative	0.09% Sodium Azide
Purity	>95%
Buffer	Borate Buffered Saline (BBS) pH 8.2 with 0.1% BSA and

Product Description


Host	Rabbit
Gene ID	970
Gene Symbol	CD70
Species	Human
Immunogen	residues 50-193 (ECD)

Product Application Details


Applications	Western Blot
Recommended Dilutions	Western Blot

Images

Western Blot: CD27 Ligand/TNFSF7/CD70 Antibody (BLR176J) [NBP3-14748] - Whole cell lysate (50 ug) from 786-O, HEK293T, HeLa, and Hs 294T cells. Rabbit anti-CD70 recombinant monoclonal antibody used at 1:1000. HRP-conjugated goat anti-rabbit IgG (A120-101P). Chemiluminescence with an exposure time of 30 seconds. Lower Panel: Rabbit anti-Actin recombinant monoclonal antibody.



Western Blot: CD27 Ligand/TNFSF7/CD70 Antibody (BLR176J) [NBP3-14748] - Whole cell lysate from HEK293F empty vector cells (10 ug) and CD70 overexpressing HEK293F cells (1 ug). Rabbit anti-CD70 recombinant monoclonal antibody used at 1:1000. HRP-conjugated goat anti-rabbit IgG (A120-101P). Chemiluminescence with an exposure time of 10 seconds. Lower Panel: Rabbit anti-Actin recombinant monoclonal antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950

Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400

Toll Free: 855.668.8722

Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15

Fax: (44) (0) 1235 533420

info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com

Technical Support: nb-technical@bio-techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Products Related to NBP3-14748

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
210-TA-005	TNF-alpha [Unconjugated]
783-CL-050/CF	CD27 Ligand/TNFSF7/CD70

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-14748

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

