

Product Datasheet

Histone Deacetylase 2/HDAC2 Antibody (S06-4E4) NBP3-14893-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-14893

Updated 4/15/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-14893



NBP3-14893-100ul

Histone Deacetylase 2/HDAC2 Antibody (S06-4E4)

Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	S06-4E4
Preservative	0.01% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 40% Glycerol, 0.05% BSA
Target Molecular Weight	60 kDa


Product Description	
Host	Rabbit
Gene Symbol	HDAC2
Species	Human, Mouse, Rat
Immunogen	Recombinant protein of human Histone Deacetylase 2/HDAC2 (Uniprot#: Q92769)

Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunofluorescence
Recommended Dilutions	Western Blot, ELISA 1:100-1:2000, Immunohistochemistry 1:500, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin, Immunofluorescence 20 ug/ml




Images


Western Blot: Histone Deacetylase 2/HDAC2 Antibody (S06-4E4) [NBP3-14893] - Western blot detection of Histone Deacetylase 2/HDAC2 in K562, C6, 3T3 cell lysates using Histone Deacetylase 2/HDAC2 Rabbit mAb (1:1000 diluted). Predicted band size: 55kDa. Observed band size: 60kDa.



Immunocytochemistry/Immunofluorescence: Histone Deacetylase 2/HDAC2 Antibody (S06-4E4) [NBP3-14893] - Immunofluorescence of Histone Deacetylase 2/HDAC2 (green) in Hela cells using Histone Deacetylase 2/HDAC2 Rabbit mAb at dilution 1/200, and DAPI (blue)



Immunohistochemistry-Paraffin: Histone Deacetylase 2/HDAC2 Antibody (S06-4E4) [NBP3-14893] - Immunohistochemistry of Histone Deacetylase 2/HDAC2 in paraffin-embedded Human colon cancer tissue using Histone Deacetylase 2/HDAC2 Rabbit mAb at dilution 1/1





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-14893-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-72399-50ug	Recombinant Human Histone Deacetylase 2/HDAC2 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-14893

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

