

# Product Datasheet

## hnRNP K/J Antibody (S07-5H0) NBP3-14943-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-14943](http://www.novusbio.com/NBP3-14943)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-14943](http://www.novusbio.com/reviews/destination/NBP3-14943)



**NBP3-14943-100ul**

hnRNP K/J Antibody (S07-5H0)

| Product Information     |                                                                                             |
|-------------------------|---------------------------------------------------------------------------------------------|
| Unit Size               | 100 ul                                                                                      |
| Concentration           | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                 | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.      |
| Clonality               | Monoclonal                                                                                  |
| Clone                   | S07-5H0                                                                                     |
| Preservative            | 0.01% Sodium Azide                                                                          |
| Isotype                 | IgG                                                                                         |
| Purity                  | Affinity purified                                                                           |
| Buffer                  | 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 40% Glycerol, 0.05% BSA                             |
| Target Molecular Weight | 62 kDa                                                                                      |


| Product Description |                                                                                                                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description         | Novus Biologicals Rabbit hnRNP K/J Antibody (S07-5H0) (NBP3-14943) is a recombinant monoclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host                | Rabbit                                                                                                                                                                                                       |
| Gene ID             | 3190                                                                                                                                                                                                         |
| Gene Symbol         | HNRNPK                                                                                                                                                                                                       |
| Species             | Human, Mouse, Rat                                                                                                                                                                                            |
| Immunogen           | A synthetic peptide of human hnRNP K/J (Uniprot#: P61978)                                                                                                                                                    |

| Product Application Details |                                                                                           |
|-----------------------------|-------------------------------------------------------------------------------------------|
| Applications                | Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry                         |
| Recommended Dilutions       | Western Blot 1/1000-1/5000, Immunohistochemistry 1/20-1/50, Immunohistochemistry-Paraffin |




## Images

Western Blot: hnRNP K/J Antibody (S07-5H0) [NBP3-14943] - Western blot detection of hnRNP K/J in K562, PC3, Hela cell lysates using hnRNP K/J Rabbit mAb (1:500 diluted). Predicted band size: 51kDa. Observed band size: 62kDa.



Immunohistochemistry-Paraffin: hnRNP K/J Antibody (S07-5H0) [NBP3-14943] - Immunohistochemistry of hnRNP K/J in paraffin-embedded Human lung cancer tissue using hnRNP K/J Rabbit mAb at dilution 1/50





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-14943-100ul**

---

|            |                                                     |
|------------|-----------------------------------------------------|
| HAF008     | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160     | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control                          |
| 210-TA-005 | TNF-alpha [Unconjugated]                            |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-14943](http://www.novusbio.com/reviews/submit/NBP3-14943)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

