

Product Datasheet

Rad51 Antibody (S04-3C6)

NBP3-14968-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-14968

Updated 2/20/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-14968



NBP3-14968-100ul

Rad51 Antibody (S04-3C6)

Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	S04-3C6
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	10mM PBS (pH 7.4), 150mM sodium chloride, 0.05% BSA, 50% glycerol
Target Molecular Weight	37 kDa


Product Description	
Host	Rabbit
Gene ID	5888
Gene Symbol	RAD51
Species	Human, Mouse, Rat
Immunogen	A synthetic peptide of human Rad51 (Uniprot#: Q06609)

Product Application Details	
Applications	Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation
Recommended Dilutions	Western Blot 1:1000-1:5000, Flow Cytometry 1:20-1:100, Immunohistochemistry 1:100-1:200, Immunocytochemistry/ Immunofluorescence 1:50-1:200, Immunoprecipitation 1:20-1:50, Immunohistochemistry-Paraffin




Images

Western Blot: Rad51 Antibody (S04-3C6) [NBP3-14968] - Western blot detection of Rad51 in A549, HL-60, U2OS cell lysates using Rad51 Rabbit mAb (1:1000 diluted). Predicted band size: 37kDa. Observed band size: 37kDa.



Immunohistochemistry-Paraffin: Rad51 Antibody (S04-3C6) [NBP3-14968] - Immunohistochemistry of Rad51 in paraffin-embedded Human tonsil using Rad51 Rabbit mAb at dilution 1/50





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-14968-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00005888-P01-10ug	Recombinant Human Rad51 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-14968

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

