

Product Datasheet

Cofilin Antibody (S06-2C0)

NBP3-15015-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-15015

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-15015



NBP3-15015-100ul

Cofilin Antibody (S06-2C0)

Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	S06-2C0
Preservative	0.01% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 40% Glycerol, 0.05% BSA
Target Molecular Weight	19 kDa


Product Description	
Description	Novus Biologicals Rabbit Cofilin Antibody (S06-2C0) (NBP3-15015) is a recombinant monoclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	1072
Gene Symbol	CFL1
Species	Human, Mouse, Rat
Immunogen	A synthetic peptide of human Cofilin (Uniprot#: P23528)

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1/1000, Immunohistochemistry 1/50, Immunocytochemistry/Immunofluorescence 1/50, Immunohistochemistry-Paraffin 5 ug/ml




Images


Western Blot: Cofilin Antibody (S06-2C0) [NBP3-15015] - Western blot detection of Cofilin in A549, MCF-7, 293T cell lysates using Cofilin Rabbit mAb (1:1000 diluted). Predicted band size: 19kDa. Observed band size: 19kDa.



Immunocytochemistry/Immunofluorescence: Cofilin Antibody (S06-2C0) [NBP3-15015] - Immunocytochemistry of Cofilin (green) in hela using Cofilin Rabbit mAb at dilution 1/50, and DAPI (blue)



Immunohistochemistry: Cofilin Antibody (S06-2C0) [NBP3-15015] - Immunohistochemical of Cofilin in Human lung cancer tissue using Cofilin antibody at dilution 1/2





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-15015-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBC1-18535	Recombinant Human Cofilin His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-15015

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

