

Product Datasheet

Smad1 [p Ser463, p Ser465] Antibody (S07-4A3) NBP3-15092-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-15092

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-15092



NBP3-15092-100ul

Smad1 [p Ser463, p Ser465] Antibody (S07-4A3)


Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	S07-4A3
Preservative	0.01% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 40% Glycerol, 0.05% BSA
Target Molecular Weight	52 kDa

Product Description	
Description	Novus Biologicals Rabbit Smad1 [p Ser463, p Ser465] Antibody (S07-4A3) (NBP3-15092) is a recombinant monoclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	4086
Gene Symbol	SMAD1
Species	Human, Mouse
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Ser463/Ser465 of human Smad1 (Uniprot#: Q15797)

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot 1/1000

Images

Western Blot: Smad1 [p Ser463, p Ser465] Antibody (S07-4A3) [NBP3-15092] - Western blot detection of Smad1 [p Ser463, p Ser465] in A549, U2OS, C6 cell lysates using Smad1 [p Ser463, p Ser465] Rabbit mAb (1:1000 diluted). Predicted band size: 52kDa. Observed band size: 52kDa.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-15092-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00004086-P01-10ug	Recombinant Human Smad1 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-15092

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

