

Product Datasheet

Hexokinase 1 Antibody (6B7Q2)

NBP3-15274-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-15274

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-15274



NBP3-15274-100ul

Hexokinase 1 Antibody (6B7Q2)


Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	6B7Q2
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol, 0.05% BSA

Product Description	
Description	Novus Biologicals Knockout (KO) Validated Rabbit Hexokinase 1 Antibody (6B7Q2) (NBP3-15274) is a recombinant monoclonal antibody validated for use in IHC, WB, ICC/IF and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	3098
Gene Symbol	HK1
Species	Human, Mouse, Rat
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human Hexokinase 1 (P19367). GLSRDFNPTATVKMLPTFVRSIPDGSEKGD FIALDLGGSSFRILRVQVNHEKNQ NVHMESEVYDTPENIVHGSGSQLFDHVAECLGDFMEKRRIKDKKLPV


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunoprecipitation, Knockout Validated
Recommended Dilutions	Western Blot 1:500 - 1:1000, Immunohistochemistry 1:50 - 1:200, Immunocytochemistry/ Immunofluorescence 1:50 - 1:200, Immunoprecipitation 1:500 - 1:1000, Immunohistochemistry-Paraffin, Knockout Validated

Images


Western Blot: Hexokinase 1 Antibody (6B7Q2) [NBP3-15274] - Western blot analysis of extracts of various cell lines, using Hexokinase 1 Rabbit mAb (NBP3-15274) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.




Immunocytochemistry/Immunofluorescence: Hexokinase 1 Antibody (6B7Q2) [NBP3-15274] - Immunofluorescence analysis of U-2 OS cells using Hexokinase 1 Rabbit mAb (NBP3-15274) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.




Immunohistochemistry-Paraffin: Hexokinase 1 Antibody (6B7Q2) [NBP3-15274] - Immunohistochemistry of paraffin-embedded mouse kidney using Hexokinase 1 Rabbit mAb (NBP3-15274) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.




Immunohistochemistry-Paraffin: Hexokinase 1 Antibody (6B7Q2) [NBP3-15274] - Immunohistochemistry of paraffin-embedded rat kidney using Hexokinase 1 Rabbit mAb (NBP3-15274) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.




Immunohistochemistry-Paraffin: Hexokinase 1 Antibody (6B7Q2) [NBP3-15274] - Immunohistochemistry of paraffin-embedded human thyroid cancer using Hexokinase 1 Rabbit mAb (NBP3-15274) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.




Western Blot: Hexokinase 1 Antibody (6B7Q2) [NBP3-15274] - Western blot analysis of extracts from normal (control) and Hexokinase 1 knockout (KO) 293T cells, using Hexokinase 1 Rabbit mAb (NBP3-15274) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.




Immunoprecipitation: Hexokinase 1 Antibody (6B7Q2) [NBP3-15274] - Immunoprecipitation analysis of 600ug extracts of Mouse brain cells using 3ug Hexokinase 1 antibody (NBP3-15274). Western blot was performed from the immunoprecipitate using Hexokinase 1 (NBP3-15274) at a dilution of 1:1000.




Immunocytochemistry/ Immunofluorescence: Hexokinase 1 Antibody (6B7Q2) [Hexokinase 1] - Confocal imaging of U-2 OS cells using [KO Validated] Hexokinase 1 Rabbit mAb followed by a further incubation with Cy3 Goat Anti-Rabbit IgG (H+L). The cells were counterstained with alpha-Tubulin Mouse mAb followed by incubation with ABflo 488-conjugated Goat Anti-Mouse IgG (H+L) Ab (Green). DAPI was used for nuclear staining (Blue). Objective: 100x.




Immunohistochemistry: Hexokinase 1 Antibody (6B7Q2) [Hexokinase 1] - Immunohistochemistry analysis of paraffin-embedded Rat kidney tissue using [KO Validated] Hexokinase 1 Rabbit mAb at dilution of 1:200 (40x lens). High pressure antigen retrieval performed with 0.01M Citrate Bufferr (pH 6.0) prior to IHC staining.



Immunocytochemistry/ Immunofluorescence: Hexokinase 1 Antibody (6B7Q2) [Hexokinase 1] - Confocal imaging of paraffin-embedded Mouse kidney tissue using [KO Validated] Hexokinase 1 Rabbit mAb followed by a further incubation with Cy3 Goat Anti-Rabbit IgG (H+L) . DAPI was used for nuclear staining (Blue). Objective: 40x. Perform high pressure antigen retrieval with 0.01M citrate buffer (pH 6.0) prior to IF staining.



Immunohistochemistry: Hexokinase 1 Antibody (6B7Q2) [Hexokinase 1] - Immunohistochemistry analysis of paraffin-embedded Mouse kidney tissue using [KO Validated] Hexokinase 1 Rabbit mAb at dilution of 1:200 (40x lens). High pressure antigen retrieval performed with 0.01M Citrate Bufferr (pH 6.0) prior to IHC staining.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-15274-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-85789PEP	Hexokinase 1 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-15274

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

