

Product Datasheet

RAP30 Antibody (3T7Q9)

NBP3-15275-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-15275

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-15275



NBP3-15275-100ul**RAP30 Antibody (3T7Q9)**


Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	3T7Q9
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol, 0.05% BSA

Product Description	
Description	Novus Biologicals Rabbit RAP30 Antibody (3T7Q9) (NBP3-15275) is a recombinant monoclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	2963
Gene Symbol	GTF2F2
Species	Human, Mouse, Rat
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RAP30 (P13984). MAERGELDLTGAKQNTGVWLVKVPKYLSQQWAKASGRGEV GKLRIAKTQGRT EVSFTLNEDLANIHDIGGKPASVSAPREHPFVLQSVGGQTLTVFTSS


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500 - 1:1000, Immunohistochemistry 1:50 - 1:200, Immunohistochemistry-Paraffin

Images


Western Blot: RAP30 Antibody (3T7Q9) [NBP3-15275] - Western blot analysis of extracts of various cell lines, using RAP30 antibody (NBP3-15275) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.




Immunohistochemistry-Paraffin: RAP30 Antibody (3T7Q9) [NBP3-15275]
 - Immunohistochemistry of paraffin-embedded rat testis using RAP30 Rabbit mAb (NBP3-15275) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.




Immunohistochemistry-Paraffin: RAP30 Antibody (3T7Q9) [NBP3-15275]
 - Immunohistochemistry of paraffin-embedded human esophageal using RAP30 Rabbit mAb (NBP3-15275) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.




Immunohistochemistry-Paraffin: RAP30 Antibody (3T7Q9) [NBP3-15275]
 - Immunohistochemistry of paraffin-embedded mouse kidney using RAP30 Rabbit mAb (NBP3-15275) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.




Immunohistochemistry: RAP30 Antibody (3T7Q9) [NBP3-15275] -
 Immunohistochemistry analysis of RAP30 in paraffin-embedded human colon tissue using RAP30 Rabbit mAb at a dilution of 1:200 (40x lens). High pressure antigen retrieval was performed with 0.01 M citrate buffer (pH 6.0) prior to IHC staining.




Immunohistochemistry: RAP30 Antibody (3T7Q9) [NBP3-15275] - Immunohistochemistry analysis of RAP30 in paraffin-embedded human lung tissue using RAP30 Rabbit mAb at a dilution of 1:200 (40x lens). High pressure antigen retrieval was performed with 0.01 M citrate buffer (pH 6.0) prior to IHC staining.




Immunohistochemistry: RAP30 Antibody (3T7Q9) [NBP3-15275] - Immunohistochemistry analysis of RAP30 in paraffin-embedded human esophagus tissue using RAP30 Rabbit mAb at a dilution of 1:200 (40x lens). High pressure antigen retrieval was performed with 0.01 M citrate buffer (pH 6.0) prior to IHC staining.



Immunohistochemistry: RAP30 Antibody (3T7Q9) [NBP3-15275] - Immunohistochemistry analysis of RAP30 in paraffin-embedded mouse testis tissue using RAP30 Rabbit mAb at a dilution of 1:200 (40x lens). High pressure antigen retrieval was performed with 0.01 M citrate buffer (pH 6.0) prior to IHC staining.



Immunohistochemistry: RAP30 Antibody (3T7Q9) [NBP3-15275] - Immunohistochemistry analysis of RAP30 in paraffin-embedded rat colon tissue using RAP30 Rabbit mAb at a dilution of 1:200 (40x lens). High pressure antigen retrieval was performed with 0.01 M citrate buffer (pH 6.0) prior to IHC staining.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-15275-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP3-17983-100ug	Recombinant Human RAP30 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-15275

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

