

Product Datasheet

Transcription factor E3 Antibody (1O9E9) NBP3-15276-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-15276

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-15276



NBP3-15276-100ul


Transcription factor E3 Antibody (1O9E9)

Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1O9E9
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol, 0.05% BSA
Product Description	
Description	Novus Biologicals Rabbit Transcription factor E3 Antibody (1O9E9) (NBP3-15276) is a recombinant monoclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	7030
Gene Symbol	TFE3
Species	Human, Mouse, Rat
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 191-440 of human Transcription factor E3 (P19532). LQQARRQQVKQYLSTTLGPKLASQALTPPPGPASQAQLPAPEAAHTTGPTGSA PNSPMALLTIGSSSEKEIDDVIDEIIISLESSYNDEMLSYLPGGTTGLQLPSTLPVS GNLLDVYSSQGVATPAITVSNSCPAELPNIKREISETAKALLKERQKKDNHNL ERRRRFNINDRIKELGTLIPKSSDPEMRWNKGTLKASVDYIRKLQKEQQRSKDL ESRQRSLEQANRSLQLRIQELELQAQIHGLP
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500 - 1:1000, Immunohistochemistry 1:50 - 1:200, Immunohistochemistry-Paraffin




Images


Western Blot: Transcription factor E3 Antibody (7V8M1) [NBP3-15276] - Western blot analysis of extracts of various cell lines, using Transcription factor E3 Rabbit mAb (NBP3-15276) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.




Immunohistochemistry-Paraffin: Transcription factor E3 Antibody (7V8M1) [NBP3-15276] - Immunohistochemistry of paraffin-embedded human esophageal cancer using Transcription factor E3 Rabbit mAb (NBP3-15276) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.



Western Blot: Transcription factor E3 Antibody (7V8M1) [NBP3-15276] - Western blot analysis of extracts of HeLa cells, using Transcription factor E3 Rabbit mAb (NBP3-15276) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry: Transcription factor E3 Antibody (1O9E9) [Transcription factor E3] - Immunohistochemistry analysis of paraffin-embedded Human acinar soft tissue sarcoma using Transcription factor E3 Rabbit mAb at dilution of 1:100 (40x lens). High pressure antigen retrieval performed with 0.01M Tris/EDTA Buffer (pH 9.0) prior to IHC staining.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-15276-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00007030-Q01-10ug	Recombinant Human Transcription factor E3 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-15276

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

