

Product Datasheet

KHSRP Antibody (4V4O1) NBP3-16746-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-16746

Updated 11/7/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-16746



NBP3-16746-100ul

KHSRP Antibody (4V4O1)

Product Information

Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	4V4O1
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS, 0.05% BSA, 50% glycerol, pH7.3

Product Description


Host	Rabbit
Gene ID	8570
Gene Symbol	KHSRP
Species	Human, Mouse, Rat
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human KHSRP (Q92945). MSDYSTGGPPPGPPPPAGGGGGAGGAGGGPPPGPPGAGDRGGGGPGGGG PGGGSAGGPSQPPGGGGPGIRKDAFADAVQRRARQIAAKIGGDAATTVNNSTP DFGFGGQKRQLEDGDQPESKKLASQGDSISSQLGPIHPPPTSMTEEYRVPD GMVGLIIGRGGEQINKIQD SGCKVQI

Product Application Details


Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500 - 1:1000, Immunohistochemistry 1:50 - 1:200, Immunoprecipitation 1:500 - 1:1000, Immunohistochemistry-Paraffin

Images


Western Blot: KHSRP Antibody (4V4O1) [NBP3-16746] - Western blot analysis of extracts of various cell lines, using KHSRP Rabbit mAb (NBP3-16746) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.




Immunohistochemistry-Paraffin: KHSRP Antibody (4V4O1) [NBP3-16746] - Immunohistochemistry of paraffin-embedded mouse kidney using KHSRP Rabbit mAb (NBP3-16746) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.




Immunohistochemistry-Paraffin: KHSRP Antibody (4V4O1) [NBP3-16746] - Immunohistochemistry of paraffin-embedded rat brain using KHSRP Rabbit mAb (NBP3-16746) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.



Immunohistochemistry-Paraffin: KHSRP Antibody (4V4O1) [NBP3-16746] - Immunohistochemistry of paraffin-embedded human esophageal cancer using KHSRP Rabbit mAb (NBP3-16746) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.



Immunoprecipitation: KHSRP Antibody (4V4O1) [NBP3-16746] - Immunoprecipitation analysis of 300ug extracts of HepG2 cells using 3ug KHSRP antibody (NBP3-16746). Western blot was performed from the immunoprecipitate using KHSRP antibody (NBP3-16746) at a dilution of 1:1000.



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-16746-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00008570-Q01-10ug	Recombinant Human KHSRP GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-16746

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

