

Product Datasheet

RAP74 Antibody NBP3-17975

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-17975

Updated 8/8/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-17975



NBP3-17975

RAP74 Antibody

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7), 20% Glycerol, 1% BSA

Product Description	
Host	Rabbit
Gene ID	2962
Gene Symbol	GTF2F1
Species	Human, Mouse
Reactivity Notes	Chimpanzee (97%), Rat (87%).
Immunogen	Recombinant protein encompassing a sequence within the center region of human TFIIIF. The exact sequence is proprietary.

Product Application Details	
Applications	Western Blot, Immunoprecipitation, Chromatin Immunoprecipitation (ChIP)
Recommended Dilutions	Western Blot 1:500-1:3000, Immunoprecipitation 1:100-1:500, Chromatin Immunoprecipitation (ChIP)



Images

Western Blot: RAP74 Antibody [NBP3-17975] - Sample (30 ug of whole cell lysate) A: 293T B: A431 C: HeLa 7. 5% SDS PAGE gel, diluted at 1:3000.




Chromatin Immunoprecipitation (ChIP): RAP74 Antibody [NBP3-17975] - Cross-linked ChIP was performed with HeLa chromatin extract and 5 ug of either control rabbit IgG or anti-TFIIF antibody. The precipitated DNA was detected by PCR with primer set targeting to GAPDH promoter or MyoD.



Western Blot: RAP74 Antibody [NBP3-17975] - Sample (30 ug of whole cell lysate) A: NIH-3T3 7. 5% SDS PAGE gel, diluted at 1:3000.



Immunoprecipitation: RAP74 Antibody [NBP3-17975] - 293T whole cell lysate/extract A : 30 ug whole cell lysate/extract of TFIIF protein expressing 293T cells B : Control with 2.5 ug of pre-immune rabbit IgG C : Immunoprecipitation of TFIIF by 2.5 ug of TFIIF antibody 7.5% SDS-PAGE The immunoprecipitated TF





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-17975

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00002962-Q01-10ug	Recombinant Human RAP74 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-17975

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

