

# Product Datasheet

## UBXD2 Antibody - BSA Free

### NBP3-18404

Unit Size: 100 ug

Store at 4C. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-18404](http://www.novusbio.com/NBP3-18404)

Updated 3/18/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-18404](http://www.novusbio.com/reviews/destination/NBP3-18404)



**NBP3-18404**

UBXD2 Antibody - BSA Free

**Product Information**

|                      |                                 |
|----------------------|---------------------------------|
| <b>Unit Size</b>     | 100 ug                          |
| <b>Concentration</b> | 1 mg/ml                         |
| <b>Storage</b>       | Store at 4C. Do not freeze.     |
| <b>Clonality</b>     | Polyclonal                      |
| <b>Preservative</b>  | 0.09% Sodium Azide              |
| <b>Isotype</b>       | IgG                             |
| <b>Purity</b>        | Affinity purified               |
| <b>Buffer</b>        | Tris-citrate/phosphate pH 7 -8. |

**Product Description**


|                    |  |
|--------------------|--|
| <b>Host</b>        | Rabbit   |
| <b>Gene ID</b>     | 23190  |
| <b>Gene Symbol</b> | UBXN4  |
| <b>Species</b>     | Human  |
| <b>Immunogen</b>   | Recognizes region between residues Between Between 458 and 508.. |

**Product Application Details**

|                              |                               |
|------------------------------|-------------------------------|
| <b>Applications</b>          | Western Blot                  |
| <b>Recommended Dilutions</b> | Western Blot 1:2000 - 1:10000 |

**Images**

Western Blot-UBXD2 Antibody-NBP3-18404- Whole cell lysate (50 ug) from HeLa, HEK293T, and Jurkat cells prepared using NETN lysis buffer. Antibody: Affinity purified rabbit anti-UBXN4 antibody used for WB at 0.1 mg/ml. Detection: Chemiluminescence with an exposure time of 30 seconds.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-18404**

---

|               |   |
|---------------|---|
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891    | Rabbit IgG Isotype Control                          |
| NBP1-92560PEP | UBXD2 Recombinant Protein Antigen                   |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-18404](http://www.novusbio.com/reviews/submit/NBP3-18404)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

