

# Product Datasheet

## DUSP12 Antibody - BSA Free NBP3-18420

Unit Size: 100 ug

Store at 4C. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-18420](http://www.novusbio.com/NBP3-18420)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-18420](http://www.novusbio.com/reviews/destination/NBP3-18420)



**NBP3-18420**

DUSP12 Antibody - BSA Free

| Product Information |                                 |
|---------------------|---------------------------------|
| Unit Size           | 100 ug                          |
| Concentration       | 1 mg/ml                         |
| Storage             | Store at 4C. Do not freeze.     |
| Clonality           | Polyclonal                      |
| Preservative        | 0.09% Sodium Azide              |
| Isotype             | IgG                             |
| Purity              | Affinity purified               |
| Buffer              | Tris-citrate/phosphate pH 7 -8. |


| Product Description |   |
|---------------------|---|
| Description         | Novus Biologicals Rabbit DUSP12 Antibody - BSA Free (NBP3-18420) is a polyclonal antibody validated for use in WB and IP. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host                | Rabbit  |
| Gene ID             | 11266   |
| Gene Symbol         | DUSP12  |
| Species             | Human   |
| Immunogen           | Recognizes region between residues Between Between 1 and 50..   |

| Product Application Details |  |
|-----------------------------|--|
| Applications                | Western Blot, Immunoprecipitation                                  |
| Recommended Dilutions       | Western Blot 1:2000 - 1:10000, Immunoprecipitation 20 ug/mg lysate |



## Images


Western Blot-DUSP12 Antibody-NBP3-18420-Whole cell lysate (50 ug) from HeLa, HEK293T, Jurkat, mouse TCMK-1, and mouse NIH 3T3 cells prepared using NETN lysis buffer. Antibody: Affinity purified rabbit anti-PITPNB antibody used for WB at 0.1 mg/ml. Detection: Chemiluminescence with an exposure time of 30 seconds.



Immunoprecipitation-DUSP12 Antibody-NBP3-18420-. Samples: Whole cell lysate (0.5 or 1.0 mg per IP reaction; 20% of IP loaded) from 293T cells prepared using NETN lysis buffer. Antibodies: Affinity purified rabbit anti-DUSP12 antibody used for IP at 2 ug per reaction. DUSP12 was also immunoprecipitated by rabbit anti-DUSP12 antibody. For blotting immunoprecipitated DUSP12, was used at 0.1 mg/ml. Detection: Chemiluminescence with an exposure time of 3 minutes.



Western Blot;DUSP12 Antibody Whole cell lysate (15  $\mu$ g) from HeLa, HEK293T, and Jurkat cells prepared using NETN lysis buffer. Antibody: Affinity purified rabbit anti-DUSP12 antibody used for WB at 0.1 mg/ml. Detection: Chemiluminescence with an exposure time of 3 minutes.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-18420**

---

|             |   |
|-------------|---|
| NBP2-33376H | Blue Marker Antibody (6F4-F6) [HRP]                 |
| HAF008      | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160      | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891  | Rabbit IgG Isotype Control                          |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-18420](http://www.novusbio.com/reviews/submit/NBP3-18420)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

