

Product Datasheet

NCBP3 Antibody - BSA Free

NBP3-18493

Unit Size: 100 ug

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-18493

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-18493



NBP3-18493

NCBP3 Antibody - BSA Free

Product Information	
Unit Size	100 ug
Concentration	1 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	Tris-citrate/phosphate pH 7 -8.


Product Description	
Description	Novus Biologicals Rabbit NCBP3 Antibody - BSA Free (NBP3-18493) is a polyclonal antibody validated for use in WB and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	55421
Gene Symbol	NCBP3
Species	Human
Immunogen	Recognizes region between residues Between between 25 and 75..

Product Application Details	
Applications	Western Blot, Immunoprecipitation
Recommended Dilutions	Western Blot 1:10000 - 1:25000, Immunoprecipitation 50-100 ug/mg lysate




Images

Western Blot;NCBP3 Antibody Whole cell lysate (50 μ g) from HeLa, LNCaP, and HEK293T cells prepared using NETN lysis buffer. Antibody: Affinity purified Rabbit anti-C17orf85 antibody used for WB at 0.04 mg/ml. Secondary: HRP-conjugated goat anti-rabbit IgG). Detection: Chemiluminescence with an exposure time of 3 seconds.



Whole cell lysate (1.0 mg per IP reaction; 20% of IP loaded) from HeLa cells prepared using NETN lysis buffer. Antibodies: Affinity purified Rabbit anti-C17orf85 antibody used for IP at 1 μ g per reaction. C17orf85 was also immunoprecipitated by rabbit anti-C17orf85 antibody BL24088. For blotting immunoprecipitated NCBP3, was used at 0.04 mg/ml. Detection: Chemiluminescence with an exposure time time of 3 seconds.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-18493

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-18493

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

