

Product Datasheet

Mouse Tenascin C ELISA Kit (Colorimetric) NBP3-18785

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-18785

Updated 1/5/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-18785



NBP3-18785**Mouse Tenascin C ELISA Kit (Colorimetric)**


| Product Information | |
|----------------------------|--|
| Unit Size | 1 Kit |
| Concentration | Concentration is not relevant for this product. Please see the protocols for proper use of this product. |
| Storage | Storage of components varies. See protocol for specific instructions. |

| Product Description | |
|--------------------------------|--|
| Gene ID | 3371 |
| Gene Symbol | TNC |
| Species | Mouse |
| Specificity/Sensitivity | This kit recognizes Mouse Tenascin C in samples. No significant cross-reactivity or interference between Mouse Tenascin C and analogues was observed. |
| Kit Components | Micro ELISA Plate (Dismountable), Reference Standard, Concentrated Biotinylated Detection Ab (100x), Concentrated HRP Conjugate (100x), Reference Standard & Sample Diluent, Biotinylated Detection Ab Diluent, HRP Conjugate Diluent, Concentrated Wash Buffer (25x), Substrate Reagent, Stop Solution, Plate Sealer, Manual, Certificate of Analysis |
| Standard Curve Range | 31.25-2000 pg/mL |
| Sensitivity | 18.75 pg/mL |
| Inter Assay | %CV < 5.09 |
| Intra Assay | %CV < 5.85 |
| Assay Type | Sandwich-ELISA |
| Spiking Recovery | 85-103% |
| Suitable Sample Type | Serum, plasma and other biological fluids |
| Sample Volume | 100 uL |

| Product Application Details | |
|------------------------------------|-----------------------|
| Applications | ELISA, Sandwich ELISA |
| Recommended Dilutions | ELISA, Sandwich ELISA |

Images

ELISA: Mouse Tenascin C ELISA Kit (Colorimetric) [NBP3-18785] - Standard Reference Curve.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to NBP3-18785

| | |
|---------------|--|
| NBP1-89639PEP | Tenascin C Recombinant Protein Antigen |
| 236-EG-200 | EGF [Unconjugated] |
| NB110-68136 | Tenascin C Antibody (4C8MS) - BSA Free |
| NBP1-91258 | Fibronectin Antibody - BSA Free |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-18785

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

