

Product Datasheet

beta-III Tubulin Antibody (S02-4G7) NBP3-19582-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-19582

Updated 4/15/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-19582



NBP3-19582-100ul


beta-III Tubulin Antibody (S02-4G7)

Product Information	
Unit Size	100 ul
Concentration	0.3 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	S02-4G7
Preservative	0.01% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.05% BSA
Target Molecular Weight	50 kDa
Product Description	
Host	Rabbit
Gene Symbol	TUBB3
Species	Human, Mouse, Rat
Immunogen	A synthetic peptide of human beta-III Tubulin (Uniprot # Q13509)
Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot, Immunohistochemistry 1:50, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin




Images


Western Blot: beta-III Tubulin Antibody (S02-4G7) [NBP3-19582] - Western blot detection of beta-III Tubulin in Rat Brain, Hela cell lysates using NBP3-19582 (1:1000 diluted). Predicted band size: 50kDa. Observed band size: 50kDa.



Immunocytochemistry/Immunofluorescence: beta-III Tubulin Antibody (S02-4G7) [NBP3-19582] - Immunofluorescence of beta-III Tubulin (green) in Hela using NBP3-19582 at dilution 1/5, and DAPI (blue)



Immunohistochemistry-Paraffin: beta-III Tubulin Antibody (S02-4G7) [NBP3-19582] - Immunohistochemistry of beta-III Tubulin in paraffin-embedded Human Cholangiocarcinoma using NBP3-19582 at dilution 1/50





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-19582-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-22901	Recombinant Human beta-III Tubulin His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-19582

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

